Non Vascular Sterile Hydro Slide Guidewire
SpiraTrax™ Non Vascular Sterile Hydro Slide Guide-
wire

Instruction for Use

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WARNING

1. The product is intended for single use only! DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed.
2. Do not use this device for any purpose other than its intended use.
3. Do not insert the device into the endoscope unless you have a clear endoscopic field of view. If you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or device.

DEVICE NAME

Non Vascular Sterile Hydro Slide Guidewire
SpiraTrax™ Non Vascular Sterile Hydro Slide Guidewire

DEVICE DESCRIPTION

【CONTENTS】


Each packaged pouch contains the following:

One (1) Guidewire

The core is made of super elastic alloy to prevent kinking. The atraumatic distal tip is made of soft polymeric material with hydrophilic coating to facilitate ease of cannulation

【CATALOGUE NUMBER】

MTN-BM-□□/□□-□-□-□-□


* If there is -W in REF, the product has a trade name "SpiraTrax™".

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Max. Diameter(mm)	Working Length(mm)	Tip Type	Shaft Type	Surface Type
≤0.45				
≤0.53	2600±130	Straight,	Normal,	Smooth,
≤0.63	4500±225	Angle	Stiff	Wave
≤0.89				

【STRUCTURE】


Figure 1 Product Drawing

No	Component	No	Component
1	Axial core	4	Hydrophilic Coating
2	Shrink Tube	5	Silicone oil
3	Tip Tube	/	/

【USER INFORMATION/TRAINING/QUALIFICATIONS】

The device is intended for use by physicians or medical personnel who have received appropriate training, in accordance with the applicable laws.

INTENDED USE/INDICATIONS

The device is indicated for use in selective cannulation of the biliary ducts including, but not limited to the common bile, cystic, right and left hepatic ducts. The device is designed to be used during endoscopic biliary procedures for catheter introduction and exchanges.

PATIENT TARGET GROUP

The device is intended for adult and adolescent populations.

CONTRAINDICATIONS

Contraindications for this device are those specific to endoscopic retrograde cholangiopancreatography (ERCP), which may include:

1. Patients who have obstruction of upper digestive tract, the device may not reach the descending duodenum;
2. Inability to pass the endoscope or tolerate the endoscopy;
3. Severe cardiopulmonary disease;

4. Severe coagulopathy;
5. Acute pancreatitis and severe biliary tract infection.
6. Severe allergic reaction to contrast medium;
7. Any other conditions the physician judges unsuitable for use.

ADVERSE EVENTS

Potential complications include:

- ◆ Pancreatitis
- ◆ Cholangitis
- ◆ Cholecystitis
- ◆ Hyperamylasemia
- ◆ Perforation
- ◆ Hemorrhage
- ◆ Sepsis
- ◆ Tissue Trauma
- ◆ Edema
- ◆ Complications which are not currently known or observed may be present.

CAUTIONS

1. Please read the instructions for use entirely before use.
2. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
3. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
4. Do not modify this device in any way, form or matter.
5. This device is not intended to be used in the presence of Flammable liquids, in an Oxygen Enriched Atmosphere or in the presence of explosive gases.

NOTES

1. A thorough understanding of the technical principles, clinical applications and risks associated with ERCP is necessary before using this device.
2. The device is intended for single use only and single person only.
3. If the device cannot be used by sudden reason, please remove the device by right way.

4. A single copy of the instructions for use is provided in a case box (2 devices). If further copies are requested, please contact Micro-Tech, we will provide free of charge.

HOW SUPPLIED

The device is supplied sterilized.

STORAGE

The device be stored in a cool dry, keep away from sunlight, and do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

COMPATIBILITY

Endoscope, Diagnostic and therapeutic devices during ERCP. Make sure that the inner diameter of catheter lumen intended for guidewire insertion is bigger than the outer diameter of guidewire, and The guide wire length is greater than the catheter length.

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1. Check the label before use, and make sure the product is used before the expiration date.
2. If the device appeared to have crack wire, broken surface, or other damage, please do not use and replace the device.
3. Make sure that the inner diameter of catheter lumen intended for guidewire insertion is bigger than the outer diameter of guidewire, and the length of guidewire is twice as long as the catheter.
4. Prior to removing the guidewire from the packaging hoop, flush with sterile water or saline. Save the hoop to store the guidewire if it will be used again during the same procedure.

Caution: Do not wipe guidewire with dry gauze.

5. The device is now ready to use.

【INSTRUCTIONS FOR USE】

This device is designed for use with an endoscope.

Caution: Do not use the guidewire if any defect is found during inspection. Please return any defective product to Micro-tech for replacement.

Caution: Do not use with metal-tip catheters. Withdrawing the guidewire through a metal tip catheter may damage surface of guidewire.

1. The device may be used in any of the following manners:
 - A. The device is placed into the duct first to gain initial access before placing a catheter over the guidewire.
 - B. A catheter is pre-loaded over the device and they are placed into the duct together.
 - C. The device is back-loaded through a pre-placed catheter.
 - D. Endoscopic devices may be exchanged over the device. The striped movement markers may be endoscopically monitored to help maintain access in the selected biliary duct. If the wire begins to move proximally or distally, the striped movement markers will provide endoscopic visual confirmation and allow repositioning of the wire.
2. Carefully advance guidewire to the desired biliary duct position while endoscopically monitoring the striped movement markers.
3. Fluoroscopically confirm the position of the distal radiopaque tip on the guidewire.
4. Complete device exchanges over the guidewire as needed by the procedure.
5. Withdraw the guidewire from the accessory channel or the device and dispose properly.

PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Manufacturer		Authorized representative in the European Community/ European Union
	Do not re-use		Do not re-sterilize
	Keep away from sunlight		Keep dry
	Use-by date		Date of manufacture
	Catalogue number		Batch code
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Sterilized using ethylene oxide		Caution
	Sterile barrier system/ sterile packaging		Open here
	Medical Device		Importer
	Unique device identifier		This device is not made with natural rubber latex

【Packaging】 Packed in flexible peel pouch

【Production Date】 See packaging

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 4 Years

WARRANTY

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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Micro-Tech (Nanjing) Co., Ltd.
MICRO-TECH Europe GmbH
Shanghai International Holding Corp. GmbH(Europe)