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WARNING:


1. The device is intended for use by physicians who have received ERCP operation training.
2. **For single use only.** Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, reprocessed or resterilized.
3. Do not use this device for any purpose other than the stated intended use.
4. Check for proper position of the balloon catheter using contrast medium injection and fluoroscopy. Balloon inflation in an improper location may lead to patient injury.
5. If resistance is met during procedure, do not advance the balloon without first determining the cause of resistance and taking remedial action.
6. Contents supplied **STERILE**. Do not use if sterile barrier is damaged. If damage is found contact the manufacturer: Micro-Tech (Nanjing) Co., Ltd.
7. The product is intended for adults.

【DEVICE NAME】

Sterile Biliary Stone Retrieval Balloon Catheter

Retrieval Balloon / short-wire compatible

【MODEL NUMBER】

(1) Product code:

SRB- Sterile Biliary Stone Retrieval Balloon Catheter

RSRB- Retrieval Balloon / short-wire compatible

(2) Three lumen catheter

T means three lumen catheter, P means three-lumen catheter with guidewire port;

(3) Balloon working diameter

(4) Working Length: 20 means 2000mm

(5) Preloaded guide wire

With-G: Preloaded guide wire; Without-G: Non Preloaded guide wire;

(6) Injection port position

With-B: Injection port is under the balloon; Without-B: Injection port is above the balloon

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Sterile Biliary Stone Retrieval Balloon Catheter

No.	REF	Work Diameter(mm) (tolerance: ± 1.5)	Working Length (mm) (tolerance: ± 100)	Matching Guide Wire Max OD(inch)	Injection Port	Channel of Endoscope (mm)
1	SRB-T-9/12-20	9/12	2000	0.035	Above	≥ 3.2
2	SRB-T-12/15-20	12/15	2000	0.035	Above	≥ 3.2
3	SRB-T-15/18-20	15/18	2000	0.035	Above	≥ 3.2
4	SRB-T-18/21-20	18/21	2000	0.035	Above	≥ 3.2
5	SRB-T-9/12/15-20	9/12/15	2000	0.035	Above	≥ 3.2
6	SRB-T-12/15/18-20	12/15/18	2000	0.035	Above	≥ 3.2
7	SRB-T-15/18/21-20	15/18/21	2000	0.035	Above	≥ 3.2

No.	REF	Work Diameter(mm) (tolerance: ± 1.5)	Working Length (mm) (tolerance: ± 100)	Matching Guide Wire Max OD(inch)	Injection Port	Channel of Endoscope (mm)
8	SRB-T-9/12-20-B	9/12	2000	0.035	Below	≥ 3.2
9	SRB-T-12/15-20-B	12/15	2000	0.035	Below	≥ 3.2
10	SRB-T-15/18-20-B	15/18	2000	0.035	Below	≥ 3.2
11	SRB-T-18/21-20-B	18/21	2000	0.035	Below	≥ 3.2
12	SRB-T-9/12/15-20-B	9/12/15	2000	0.035	Below	≥ 3.2
13	SRB-T-12/15/18-20-B	12/15/18	2000	0.035	Below	≥ 3.2
14	SRB-T-15/18/21-20-B	15/18/21	2000	0.035	Below	≥ 3.2

Retrieval Balloon / short-wire compatible


No.	REF	Work Diameter(mm) (tolerance: ± 1.5)	Working Length (mm) (tolerance: ± 100)	Matching Guide Wire Max OD(inch)	Injection Port	Channel of Endoscope (mm)
1	RSRB-T-9/12-20	9/12	2000	0.035	Above	≥ 3.2
2	RSRB-T-12/15-20	12/15	2000	0.035	Above	≥ 3.2
3	RSRB-T-15/18-20	15/18	2000	0.035	Above	≥ 3.2
4	RSRB-T-18/21-20	18/21	2000	0.035	Above	≥ 3.2
5	RSRB-T-9/12/15-20	9/12/15	2000	0.035	Above	≥ 3.2
6	RSRB-T-12/15/18-20	12/15/18	2000	0.035	Above	≥ 3.2
7	RSRB-T-15/18/21-20	15/18/21	2000	0.035	Above	≥ 3.2
8	RSRB-T-9/12-20-B	9/12	2000	0.035	Below	≥ 3.2
9	RSRB-T-12/15-20-B	12/15	2000	0.035	Below	≥ 3.2
10	RSRB-T-15/18-20-B	15/18	2000	0.035	Below	≥ 3.2
11	RSRB-T-18/21-20-B	18/21	2000	0.035	Below	≥ 3.2
12	RSRB-T-9/12/15-20-B	9/12/15	2000	0.035	Below	≥ 3.2
13	RSRB-T-12/15/18-20-B	12/15/18	2000	0.035	Below	≥ 3.2
14	RSRB-T-15/18/21-20-B	15/18/21	2000	0.035	Below	≥ 3.2
15	RSRB-P-9/12-20	9/12	2000	0.035	Above	≥ 3.2
16	RSRB-P-12/15-20	12/15	2000	0.035	Above	≥ 3.2
17	RSRB-P-15/18-20	15/18	2000	0.035	Above	≥ 3.2
18	RSRB-P-18/21-20	18/21	2000	0.035	Above	≥ 3.2
19	RSRB-P-9/12/15-20	9/12/15	2000	0.035	Above	≥ 3.2
20	RSRB-P-12/15/18-20	12/15/18	2000	0.035	Above	≥ 3.2
21	RSRB-P-15/18/21-20	15/18/21	2000	0.035	Above	≥ 3.2
22	RSRB-P-9/12-20-G	9/12	2000	0.035	Below	≥ 3.2
23	RSRB-P-12/15-20-G	12/15	2000	0.035	Below	≥ 3.2
24	RSRB-P-15/18-20-G	15/18	2000	0.035	Below	≥ 3.2
25	RSRB-P-18/21-20-G	18/21	2000	0.035	Below	≥ 3.2
26	RSRB-P-9/12/15-20-G	9/12/15	2000	0.035	Below	≥ 3.2
27	RSRB-P-12/15/18-20-G	12/15/18	2000	0.035	Below	≥ 3.2
28	RSRB-P-15/18/21-20-G	15/18/21	2000	0.035	Below	≥ 3.2
29	RSRB-P-9/12-20-B	9/12	2000	0.035	Above	≥ 3.2
30	RSRB-P-12/15-20-B	12/15	2000	0.035	Above	≥ 3.2
31	RSRB-P-15/18-20-B	15/18	2000	0.035	Above	≥ 3.2
32	RSRB-P-18/21-20-B	18/21	2000	0.035	Above	≥ 3.2

No.	REF	Work Diameter(mm) (tolerance: ± 1.5)	Working Length (mm) (tolerance: ± 100)	Matching Guide Wire Max OD(inch)	Injection Port	Channel of Endoscope (mm)
33	RSRB-P-9/12/15-20-B	9/12/15	2000	0.035	Above	≥ 3.2
34	RSRB-P-12/15/18-20-B	12/15/18	2000	0.035	Above	≥ 3.2
35	RSRB-P-15/18/21-20-B	15/18/21	2000	0.035	Above	≥ 3.2
36	RSRB-P-9/12-20-G-B	9/12	2000	0.035	Below	≥ 3.2
37	RSRB-P-12/15-20-G-B	12/15	2000	0.035	Below	≥ 3.2
38	RSRB-P-12/15-20-G-B	15/18	2000	0.035	Below	≥ 3.2
39	RSRB-P-18/21-20-G-B	18/21	2000	0.035	Below	≥ 3.2
40	RSRB-P-9/12/15-20-G-B	9/12/15	2000	0.035	Below	≥ 3.2
41	RSRB-P-12/15/18-20-G-B	12/15/18	2000	0.035	Below	≥ 3.2
42	RSRB-P-15/18/21-20-G-B	15/18/21	2000	0.035	Below	≥ 3.2

【STRUCTURE】

Sterile Biliary Stone Retrieval Balloon Catheter

The Sterile Biliary Stone Retrieval Balloon Catheter consists of Balloon, Catheter and Handle. The Handle includes the Balloon lumen, Guide-wire lumen and Injection lumen (Fig.1).


1-Balloon; 2-Catheter; 3-Handle; 4-Balloon lumen; 5-Guide-Wire lumen; 6-Injection lumen;

Fig.1 Sterile Biliary Stone Retrieval Balloon Catheter

Note: This is a sketch map of Sterile Biliary Stone Retrieval Balloon Catheter, it is for reference only, as the actual device may vary.

Retrieval Balloon / short-wire compatible (no guidewire port)

Retrieval Balloon / short-wire compatible (no guidewire port) consists of Balloon, Catheter and Handle. The Handle includes Balloon lumen and Injection lumen. There is a Rapid Exchange Ramp on the Catheter. The guidewire channel of the catheter accommodates a 0.035 in (0.89mm) guidewire from the tip of the Catheter to the Rapid Exchange Ramp. (Fig.2).


1-Balloon; 2-Catheter; 3-Handle; 4-Balloon lumen; 5- Injection lumen;

Fig.2 Retrieval Balloon / short-wire compatible (no guidewire port)

Note: This is a sketch map of Retrieval Balloon / short-wire compatible (no guidewire port), it is for reference only, as the actual device may vary.

Retrieval Balloon / short-wire compatible (with guidewire port)

Retrieval Balloon / short-wire compatible (with guidewire port) consists of Balloon, Catheter and Handle. The Handle includes Balloon lumen and Injection lumen. There is a Guidewire Port on the Catheter. The guidewire channel of the Catheter accommodates a 0.035 in (0.89mm) guidewire from the tip of the Catheter to Guidewire Port. (Fig.3).

1-Balloon; 2-Catheter; 3-Handle; 4-Balloon lumen; 5- Injection lumen; 6- Guidewire Port;

7- Guidewire (only for preloaded guidewire)

Fig.3 Retrieval Balloon / short-wire compatible (with guidewire port)

【USER INFORMATION/TRAINING/QUALIFICATIONS】

According to law, this instrument is intended for use by physicians who have received ERCP operation training.

【RADIOPACITY INFORMATION】

The Device contains Radiopaque Marker Band and is visible on X-ray.

【INTENDED USE / INDICATIONS】

The Sterile Biliary Stone Retrieval Balloon Catheter is used for the endoscopic removal of biliary stones or to facilitate the injection of contrast medium while occluding the duct with the balloon.

【CONTRAINDICATIONS】

Contraindications for this device are those for endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic sphincterotomy (ES), include:

- Patients who have obstruction of upper digestive tract, the device may not reach the descending duodenum;
- Inability to pass the endoscope or tolerate the endoscopy;
- Severe cardiopulmonary disease;
- Severe coagulopathy;
- Acute pancreatitis and severe biliary tract infection.
- Severe allergic reaction to contrast medium;
- Any other conditions the physician judges unsuitable for use.

【POTENTIAL COMPLICATIONS】

- Pancreatitis
- Cholangitis
- Hyperamylasemia
- Perforation
- Hemorrhage
- Sepsis
- Stone impaction
- Complications which are not currently known or observed may be present.

【CLINICAL BENEFIT】

The Stone Retrieval Balloon is used for the endoscopic removal of biliary stones or to facilitate the injection of contrast medium while occluding the duct with the balloon, which can relieve biliary obstruction.

【PERFORMANCE CHARACTERISTICS】

The device has an injection lumen, which can inject contrast medium.

The device can enter the human bile duct along the guidewire through the duodenoscope.

The balloon can be inflated to two or three different diameters.

【CAUTIONS】

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
3. Please read this instruction manual entirely before use.
4. Refer to package label for minimum channel size required for this device.

This product contains **NATURAL RUBBER LATEX**, which may cause allergic reactions.

【NOTES】

A single copy of the instructions for use is provided in a case box. If further copies are requested, please contact Micro-Tech, we will provide free of charge.

【STORAGE】

1. Storage conditions: in a room with a temperature range of 10-25°C and a humidity range of $\leq 80\%$.
2. Rotate inventory so that products are used prior to the sterilization expiration date on the package label.
3. Store the device in a dark place, and do not expose the package to organic solvents, ionizing radiation or ultraviolet light.

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1. Before using, inspect the pouch for any breach of the package. Do not use if package is damaged.
2. After verifying the sterilization expiry date, open the package carefully, take out the product.
3. Pull the Sterile Biliary Stone Retrieval Balloon Catheter out of the protecting tube, and then remove the protective caps on the Luer tapers (if applicable).
4. If this device shows any sign of damage and suspicious problems, do not use it. Do not attempt to repair nonfunctional or damaged devices.
5. In preparation for use, after pulling the syringe to the maximum mark, attach a stopcock to the syringe and inflate the balloon with air. Turn off the stopcock and confirm the balloon remains inflated. If leakage is observed, do not use the device, and contact Micro-Tech for replacement.

Note: Stopcock is in open position, allowing access to balloon, when two arms are aligned with catheter and syringe. To maintain balloon inflation, turn stopcock arm 90°.

【INSTRUCTIONS FOR USE】

Instruction for use of Sterile Biliary Stone Retrieval Balloon Catheter

Note: Use guide wire with 0.035 inch (0.89mm).

1. Carefully introduce the deflated balloon into the accessory channel of the duodenoscope over the previously placed guidewire. Advance the device in small increments until it appears in endoscopic view.
2. Using fluoroscopic imaging, position the deflated balloon in the duct above the stone, noting the location of the fluoroscopic marker bands on the balloon catheter.

Note: If removing more than one stone, it is recommended that they are removed one at a time.

3. Verify desired position and inflate the balloon using the selected syringe, with air only. Inflate the balloon to occlude the duct, as visualized under fluoroscopic imaging and lock the stopcock to maintain the balloon size.

Note: For the 9mm size of the 9-12mm balloon, inflate to the 12mm marker first, then to the 9mm marker on the syringe.

4. Contrast may now be injected through the injection port with wire guide in place if desired.

5. Using fluoroscopic visualization and keeping endoscope elevator open, gently withdraw inflated balloon toward papilla.

Note:

- (1) Do not exert excessive pressure on ampulla while extracting stones. If stone does not pass easily, reassess need for sphincterotomy.
- (2) In the stone removal process, when the large-diameter balloon runs from the thick bile duct to the fine bile duct, there will be an increase in drag resistance. This moment, the air volume in the balloon can be adjusted by controlling the two-way valve to reduce the diameter of the balloon to adapt to the different diameter segments of the bile duct.

6. When the balloon can be endoscopically visualized within the duodenum, open the stopcock, and allow the balloon to deflate.

7. Repeat the above procedure, until all the stones are cleared.

Note: If desired the previously placed guidewire may be left in position while the balloon catheter is removed, to facilitate the introduction of additional wire guided devices. Remove the balloon using standard device exchange techniques.

Instruction for use of Retrieval Balloon / short-wire compatible

Note: The device may be used with a Guidewire Locking Device, Use guide wire with 0.035 inch (0.89mm).

If using the device with prepositioned guidewire.

1. Insert guidewire from the distal end of the device. Slide the device over the guidewire, and separate the guide wire from the slot in the side wall of the catheter 15-50 cm from the catheter tip (With guidewire port specifications), or the guide wire is threaded through the rapid exchange ramp 40-50cm from the tip of the catheter (Non-guidewire port specifications). Holding the guide wire and slide the device over the guidewire until the tip of the catheter reach the Guidewire Locking Device
2. Insert the distal 15-50cm of the catheter into the scope. Lock the guide wire on the Guidewire locking Device and continue advancing the device though the scope until the distal portion of the catheter exits from the scope
3. Position deflated balloon above the stone(s) be retrieved from the appropriate duct

Note:

- (1) Contrast may be injected through injection port; radiopaque markers located just proximal to the balloon will be aid in positioning the catheter while using fluoroscopy
- (2) If more than one stone is to be removed, extract one stone at a time.
4. Inflate balloon with air, then lock the stopcock.

Note: for multiple sizing balloons, inflate as indicated on syringe

5. Using fluoroscopic visualization and keeping endoscope elevator open, gently withdraw inflated balloon toward papilla and slowly retrieve the stone.

Note:

- (1) Do not exert excessive pressure on ampulla while extracting stones. If stone does not pass easily, reassess need for sphincterotomy. In the process of taking the stone, it is possible to inject the contrast agent into the bile duct to check the stone has been removed completely or not under the X-ray

- (2) In the stone removal process, when the large-diameter balloon runs from the thick bile duct to the fine bile duct, there will be an increase in drag resistance. This moment, the air volume in the balloon can be adjusted by controlling the two-way valve to reduce the diameter of the balloon to adapt to the different diameter segments of the bile duct.
6. Deflate the balloon, repeat step 3, step 4 and step 5, remove the stones one by one until duct is clear.
7. Locking the guidewire on the Guidewire Locking Device.
8. Withdraw the device over the guidewire from the scope until end of the slot emerges in the device entrance of Guidewire Locking Device.
9. Unlock the guidewire and perform the standard exchange procedure until the guidewire is exposed, re-lock the guide wire on the locking device and pull the device off the guidewire.

If using the device with preloaded guide wire

10. Prior to use, adjust the guidewire so that the desired length of wire is extended beyond the device.
11. Carefully introduce the guidewire and deflated the device into the Guidewire Locking Device channel and enter the endoscopic cavity. Insert the instrument at a distance of 2-3cm each time until it appears in endoscopic view.
12. Follow steps 3-9 above

【PRODUCT DISPOSAL】

This product may be a potential biohazard after use. Handle and dispose of the product in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

【POST-PROCEDURE】

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

【HOW SUPPLIED】

The device is supplied STERILE.

COMPATIBILITY

【Applicable Endoscopes】

Endoscope which is legal listing in European Union is recommended, working channel of

Endoscope $\geq 3.2\text{mm}$

【Compatible Guidewire Diameter】

Guidewire Diameter $\leq 0.035''$

【Guidewire Locking Device】

Only applicable to Retrieval Balloon / short-wire compatible

Guidewire Locking Device for use with endoscopes with a working diameter of 3.2mm or more

SYMBOL INDICATIONS

	Do not re-use		Do not resterilize
	Date of manufacture		Manufacturer
	Use-by date		Caution
	Catalogue number		Batch code
	Consult instructions for use or consult electronic instructions for use		Sterilized using ethylene oxide
	Authorized representative in the European Community		Keep away from sunlight
	Keep dry		Do not use if package is damaged
	Medical device		Sterile barrier system / Sterile packaging
	Importer		Contains or presence of natural rubber latex
	Unique device identifier		Open here

【Packaging】 Packed in flexible peel pouch

【Production Date】 See packaging

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 Eighteen months

WARRANTY

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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MICRO-TECH Europe GmbH
Shanghai International Holding Corp. GmbH(Europe)