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WARNINGS

1. **The product is intended for single use only!** DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization, or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Re-use, re-sterilization, or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized, or reprocessed.
2. Do not use this instrument for any purpose other than its intended use.
3. Do not use this device on a patient with an implanted pacemaker, as it may cause malfunction or failure of the pacemaker, resulting in serious harm to the patient. Before proceeding, confirm with a cardiologist or the manufacturer of the pacemaker that it is safe to use this device in the presence of a pacemaker.
4. When using the instrument in the vicinity of the heart, be sure to use it with the minimum necessary output.
5. When using an electrocardiograph or other physiological monitoring equipment simultaneously with the instrument, monitoring electrodes should be placed as far away as possible from the electrodes used with the electrode unit. Needle monitoring electrodes should not be used. Physiological monitoring equipment incorporating high-frequency current limiting devices is recommended.
6. When endoscopes are used with the device, the patient leakage currents may be additive.
7. Prepare monitoring and rescue equipment to control unpredictable risks.
8. Ensure that there are no sharp edges on the device before inserting to the patients.
9. Do not modify this equipment.
10. Avoid use of highly explosive gases near the high frequency endoscopic accessories area, as it is an explosion hazard.
11. When using this device with electricity, do not touch the chute of handle (Figure A.).
12. When using this device with electricity, do not allow external conductive objects to contact the handle (Figure B.).

Figure A

Figure B

13. The patient leads should be positioned in such a way that contact with the patient or other leads is avoided.

Temporarily unused active electrodes should be stored in a location that is isolated from the patient.

14. Patients may feel a neuromuscular stimulus when there is a spark discharge to the metallic object. To

prevent this, minimize the discharge by selecting a lower output setting or by activating output when the forceps are in contact with the tissue to be cauterized or coagulated, or when obtaining hemostasis.

DEVICE NAME

【COMMON NAME】

Single-Use Coagulation Forceps

【TRADE NAME】

Ensure™ Single-Use Coagulation Forceps


DEVICE DESCRIPTION

【MODEL NUMBER/SPECIFICATION】

Unit: mm

REF	Working length(L)	Jaws type	Maximum diameter(D)	Applicable channel diameter
CF165-A	1650	Type A	≤ 2.7	≥ 2.8
CF165-B	1650	Type B	≤ 2.7	≥ 2.8
CF165-C	1650	Type C	≤ 2.7	≥ 2.8
CF230-A	2300	Type A	≤ 2.7	≥ 2.8
CF230-B	2300	Type B	≤ 2.7	≥ 2.8
CF230-C	2300	Type C	≤ 2.7	≥ 2.8

【STRUCTURE】

Structural drawing of Ensure™ Single-Use Coagulation Forceps

1. Jaws 2. Spring Sheath 3. HF Plug 4. Handle


Figure 1

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by trained health care professional only. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

INTENDED USE

These instruments have been designed to be used with endoscopes to cauterize and coagulate or to perform hemostasis using high-frequency current within the digestive tract.

INDICATIONS

Bleeding during digestive endoscopic surgery.

CONTRAINDICATIONS

Contraindications include, but are not limited to:

1. Severe cardiopulmonary insufficiency.
2. The patient's inability to withstand endoscopy procedures.
3. Severe coagulopathy.

ADVERSE EVENTS

Potential complications include, but are not limited to, rebleeding, perforation, transmural burns, thermal injury to the patient, and explosion.

CAUTIONS

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before use.
2. Inform the patient of the details of the operation and all the potential risks and complications, which may lead to injury, illness, or death, and ensure that the patient has expressed his/her acceptance.

3. Please read the instructions for use in their entirety before use.
4. The product is intended for adult populations.

HOW SUPPLIED

The Ensure™ Single-Use Coagulation Forceps is supplied STERILE.

STORAGE

The product should be stored in a clean, well-ventilated, and non-corrosive gas environment.

Do not expose the package to organic solvent, ionizing radiation, or ultraviolet radiation.

The product's shelf life is 2 years.

COMPATIBILITY

【Applicable Electrosurgical Unit】

High-frequency electrosurgical unit which is CE Marked in EUROPE is recommended, such as Erbe VIO® 200 D and its matched output line 20192-117, ConMed Electrosurgery 60-8200-230.

【Rated High-Frequency Voltage】 COAG: 2300 Vp (4600 Vp-p)

DO NOT use higher repeated peak voltage than the Coag -2300 Vp (4600 Vp-p)

【Applicable Endoscopes】 Endoscope which is CE Marked in EUROPE is recommended, such as an endoscope from Olympus.

OPERATION ENVIRONMENT

Operation environment	Environment temperature	10-40 °C (50-104 °F)
	Relative humidity	30%-85%
	Barometric pressure	700-1060hPa (0.7-1.1Kgf/cm ³) (10.2-15.4psia)

TRANSPORTATION AND STORAGE ENVIRONMENT

Transportation and storage environment	Temperature limitation	-40-+70 °C (-40-+158 °F)
	Humidity	10%-95%
	Atmospheric pressure limitation	700-1060hPa (0.7-1.1Kgf/cm ³) (10.2-15.4psia)

DIRECTIONS FOR USE

【PREPARATIVE CHECK AND PREPARTION】

1. Choose the appropriate specification by referring to the compatible working channel in the package label.
2. Contents supplied STERILE.
3. Inspect the package for damage before use. Do not use if package is damaged. Open the package carefully after verifying the period of validity.
4. Ensure that there are no sharp edges before inserting the device into the patient.
5. If this device shows any signs of damage, do not use. Do not attempt to repair a nonfunctional or damaged device.
6. Before using this device, follow recommendations provided by the electrosurgical unit manufacturer to ensure patient safety through proper placement and utilization of the patient return electrode. Ensure a proper path from the patient return electrode to electrosurgical unit is maintained throughout the procedure.

【INSTRUCTIONS FOR USE】

1. Attach the neutral electrode to the patient.
2. Endoscopically visualize the area to be cauterized or coagulated or to perform hemostasis.
3. With jaws closed, insert forceps into the accessory channel of the endoscope.
4. Advance forceps in 1cm-2 cm increments until it is visualized exiting the endoscope.

Note: Keep the end of the forceps that is extending from the accessory channel straight at all times. Allowing forceps to hang from the accessory channel may cause damage to the device.
5. Advance forceps to desired hemostasis site, then open the jaws.
6. Turn the handle to orient the jaws in the optimum direction (see Figure 2).
7. Plug the active cord plug into the plug until it clicks.
8. Following electrosurgical unit manufacturer's instructions for settings, verify desired settings and activate electrosurgical unit.
9. Using slight pressure on handle, close jaws around the tissue. After closing the jaws into the tissue and coagulating it, open the jaws again to ensure the effectiveness of coagulation, and retract the jaws from the tissue.

Caution: When applying current, ensure metal tip of jaws does not come in contact with endoscope. Contact of forceps tip with endoscope may result in grounding or injury to patient and/or operator, as well as damage to endoscope and/or forceps.
10. Advance into tissue to cauterize and coagulate or to perform hemostasis.

11. Upon completion of the procedure, turn the electrosurgical unit off and disconnect the active cord from the handle.
12. Continue to apply slight pressure on the handle and withdraw the forceps from the channel. While withdrawing the forceps from the endoscope, wipe excess secretions from the Spring Sheath.


Figure 2

【NOTES】

1. If the forceps cannot be rotated smoothly, push and pull the insertion portion little by little while turning the control section. This will make it possible to rotate the forceps smoothly.
2. If the forceps cannot be rotated smoothly, set the insertion portion protruding from the biopsy valve as straight as possible and turn the control section to retry.

EMC CONDITIONAL

Guidance and manufacturer's declaration - electromagnetic emissions - for all ME EQUIPMENT and ME SYSTEM.

Table 1 - Emission		
Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class A	Professional healthcare facility environment
Harmonic distortion	IEC 61000-3-2 Class A	
Voltage fluctuations and flicker	IEC 61000-3-3 Complies	

Guidance and manufacturer's declaration – electromagnetic immunity –
for all ME EQUIPMENT and ME SYSTEM.

Table 2 - Enclosure Port		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	3V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF	IEC 61000-4-3	Refer to Table 3

wireless communications equipment		
Electrical fast transients/bursts	IEC 61000-4-4	$\pm 2\text{kV}$, 100kHz (AC power port) $\pm 1\text{kV}$, 100kHz (signal input/output parts port)
Surges	IEC 61000-4-5	$\pm 0.5\text{kV}$, $\pm 1\text{kV}$ (line to line); $\pm 0.5\text{kV}$, $\pm 1\text{kV}$, $\pm 2\text{kV}$ (line to ground)
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz; 6V in ISM band between 0.15MHz and 80MHz
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz
Voltage dips	IEC 61000-4-11	0% U_T ; 0.5cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°; 0% U_T ; 1cycle at 0° 70% U_T ; 25/30 cycles at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T ; 250/300 cycles

Guidance and manufacturer's declaration – electromagnetic immunity –
for ME EQUIPMENT and ME SYSTEM that are not LIFE-SUPPORTING

Table 3 – Proximity fields from RF wireless communications equipment		
Test frequency (MHz)	Band (MHz)	Immunity test levels
		Professional healthcare facility environment
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, $\pm 5\text{kHz}$ deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810	800-960	Pulse modulation 18Hz, 28V/m
870		
930		
1720	1700-1990	Pulse modulation 217Hz, 28V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217Hz, 28V/m
5240	5100-5800	Pulse modulation 217Hz, 9V/m
5500		
5785		


PRODUCT DISPOSAL


After use, this product may be a potential biohazard. Handle and dispose of the device in accordance with hospital, local and administrative laws and regulations.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Diameter		Working length
	Compatible working channel		Batch code
	Catalogue Number		Authorised representative in the European Community/European Union
	Date of Manufacture		Do not re-use
	Manufacturer		Caution
	Sterilization using ethylene oxide		Keep away from sunlight
	Use-by date		Keep dry
	Do not resterilize		Type BF applied part
	Do not use if package is damaged and consult instruction for use		CE sign with the number of the notified body
	Consult instructions for use or consult electronic instructions for use		Medical device
	This device does not made with natural rubber latex.		Sterile barrier system/sterile packaging
	Open here		Temperature limitation

	Humidity limitation		Atmospheric pressure limitation
	Importer		Unique device identifier

【Packaging】 Flexible peel pouch

【Production Date】 See packing

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 Two years

Warranty

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp. GmbH (Europe)
Lotus Global Co., Ltd.