Sterile Hot Snare
Instructions for Use

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
WARNING

1. **For single use only.** Do not reuse, reprocess or re-sterilize. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to devices reused, reprocessed or re-sterilized.
2. The device, when applied to a patient with a pacemaker implanted, may cause malfunction or failure of the pacemaker, seriously affecting the patient. Consulting cardiologist or the manufacturer of the pacemaker is recommended before device application.
3. When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.
4. When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device, any monitoring electrodes should be placed as far away as possible from the electrodes used with the electrodes unit. Needle monitoring electrodes should not be used. Physiological monitoring equipment incorporating high-frequency current limiting devices is recommended.
5. When the endoscopes are used with the device, the leakage currents may be relatively increased. Ensure a proper path from patient return electrode to electrosurgical unit is connected and maintained throughout the procedure.
6. Monitoring equipment and rescue equipment are recommended for any unpredictable risks.
7. Before use, please check the insertion patients body parts, ensure that no sharp edge.
8. No modification of this equipment is allowed.
9. This device in the process of the electricity, hands cannot touch the chute of handle, like picture a.
10. This device in the process of the electricity, Shall not make external conductive objects in contact with the handle, like picture b.

picture a

picture b

11. Do not insert the device into the endoscope unless a clear endoscopic field of view. If you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. It may also damage the endoscope and/or device.
12. Since supporting the use of Olympus endoscopic system, Olympus endoscopic system's light source may affect Sterile Hot Snare's temperature to exceed 41°C, to minimize the risk of injury, please ensure the correct use of it according to Olympus endoscopic system and the manufacturer's instructions.

DEVICE NAME

Sterile Hot Snare

DEVICE DESCRIPTION

【CONTENTS】

Each snare is made of a specially stranded wire, with en-during radial force, high dimensional stability and precise cutting properties. The 3-ring handle offers a secure hold and transmits every movement reliably to the snare

【MODEL NUMBER】

HS XX – X XX XX XXX
 (1) (2) (3) (4) (5)

(1) Loop wire rope structure

01: Snare loop applies 1*7 rope

02: Snare loop applies 3*3 rope

(2) The shape of loop

0: Oval 1: Hexagonal 2: Crescent 3: Combinational 5: Diamond

(3) The width of loop

10:10mm	15:15mm	20:20mm	24:24mm
25:25mm	30:30mm	36:36mm	

(4) The diameter of outer tube

23: 2.3mm

(5) Working length of Sterile Hot Snare

160: 1600mm	180: 1800mm
230: 2300mm	250: 2500mm

For example, HS01-02423230 Snare loop applies 1*7 rope, oval snare, width of loop is 24mm, outer tube diameter 2.3mm, working length 2300mm device.

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
The specifications of the device are represented by loop width, diameter of outer tube and working length, as indicated in Table 1 .

Table 1 Specification and parameter

NO.	Snare Loop Shape	Loop width (mm)	Diameter of outer tube (mm)	Working length (mm)
1	Oval	10, 15, 20 24, 30, 36	ϕ 2.3	1600, 1800 2300, 2500
2	Diamond	10, 15 20, 25	ϕ 2.3	1800, 2300
3	Hexagonal	15, 20, 24	ϕ 2.3	1800, 2300
4	Combinational	30	ϕ 2.3	1800, 2300
5	Crescent	24	ϕ 2.3	1800, 2300

【STRUCTURE】

Sterile Hot Snare mainly includes a loop assembly, outer tube assembly, HF conductive plug, handle assembly . (Figure.1)

1. Loop assembly 2. Outer tube assembly 3. HF Conductive plug 4. Handle assembly

Figure 1 Schematic diagram of a Sterile Hot Snare

Note: This is the sketch map of one kind of loop, the other loop shapes are not shown here

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by physicians or medical personnel who have received appropriate training, in accordance with the applicable laws. As an accessory of digestive endoscopy, the device shall be used by professionals familiar with the operation technique of digestive endoscopy.

INTENDED USE/ INDICATIONS

This device is used in the endoscope to remove polyps or tissues in the gastrointestinal tract through high-frequency current.

CONTRAINDICATIONS

Contraindications for these devices are those specific to endoscopic polypectomy and tissue resection, which may include:

1. Inability to pass the endoscope or tolerate the endoscopy;
2. Severe cardiopulmonary disease;
3. Severe coagulopathy;
4. Polyp cancer, who has infiltrated and worsened.
5. Patients with implanted pacemaker or implanted with a metal valve / valves if cauterization is used (It is classified as relative contraindication);
6. Any other conditions the physician judges unsuitable for use.

POSSIBLE COMPLICATIONS

Possible complications include:

- 1.Immediate or delayed hemorrhage;
- 2.Perforation;
- 3.Infection;
- 4.Thermal injury;
- 5.Arrhythmia;
- 6.Complications which are not currently known or observed may be present.

CAUTIONS

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before use.
2. The patient should be informed and express acceptance of potential risks and complications of the procedure, which may lead to injury, illness, or death of the patients.
3. Please read this instruction for use in their entirety before use.

4. Serious incident relation to the device should be reported to the manufacturer and the competent authority of the Member State.
5. The device is intended for adult and adolescent populations.

NOTE

1. The device is intended for single use only. Do not re-use, or cross use it.
2. If polyps cannot be removed, and don't force operation, please change snare or use other methods.
3. For multiple polyps patients, after repeated procedures, please observe whether loop is distorted or not, if distortion, please use another one.
4. A single copy of the instructions for use is provided in a case box (10 devices). If further copies are requested, please contact Micro-Tech, we will provide free of charge.

HOW SUPPLIED

This device is for one-time use, and the device is supplied sterile.

ENVIRONMENT

Operation environment

Temperature limitation: 10 °C ~40 °C

Humidity limitation: 30%~85%

Atmospheric pressure limitation: 700hpa ~1060hpa

Transportation and storage environment

Temperature limitation: -40 °C ~60 °C

Humidity limitation: 30%~85%

Atmospheric pressure limitation: 700hpa ~1060hpa

STORAGE

1. Do not expose to organic solvents, ionizing radiation, or ultraviolet light.
2. Rotate inventory so that devices are used prior to the sterilization expiration date on the package label.
3. Store the device in a cool, dark and dry place. The sterilization is valid for 5 years.

COMPATIBILITY

【Applicable Electrosurgical Unit】

High-frequency electrosurgical units that are legal listed in EUROPE is recommended, such as CONMED Electrosurgery 60-8200-230 and its matched 474-L active cord.

【Rated High-Frequency Voltage】

CUT:1200Vp (2400Vp-p)

DO NOT use higher repeated peak voltage than the Cut- 1200Vp (2400Vp-p) .

【Applicable Endoscopes】


Endoscope that are legal listed in EUROPE are recommended such as Olympus.

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1. Inspect the package before use. The device cannot be used if the package is damaged.
2. After opening the package, operate the slider back and forth to observe whether the Sterile Hot Snare loop extends and retracts smoothly, whether there are distortion and deformation, unraveled wires, disconnections, sharp protrusions in the Sterile Hot Snare loop, and whether there are scratches or broken areas, etc. on the surface of the tube; do not use the device if any above is found.
3. Before using the device, inspect the endoscope and its accessories to be inserted into human body, to ensure there are no rough surfaces, sharp edges or protrusions which may cause safety hazard.
4. Before operation, make sure there is no inflammable gas in human body; otherwise, there is a risk of explosion. Inflammable gas must be removed before operation and / or use CO₂ to blow and spray the affected body cavity for the removal of the inflammable gas. During operation, make sure that the anesthetic, skin detergent and disinfectant are non-flammable.
5. Attach the negative plate of the high-frequency electrosurgical unit on the patient's vascular / muscle-rich part (as shown in the following figures).

【INSTRUCTIONS FOR USE】

1. The Sterile Hot Snare is capable to an endoscope with a channel of a minimum channel diameter of 2.8 mm
2. Operate the snare handle to pull the loop back into the outer tube (as shown in figure 2). The snare enters the digestive tract through the endoscopic channel hole. Under the endoscopic monitoring, the outer tube should be guaranteed to leak out of the distal end of the endoscopic channel (as shown in figure 3) :


Figure. 2


Figure. 3

3. Connect the plug of Sterile Hot Snare to the socket of high-frequency electrosurgical unit. Pay attention to the correct connection (as shown in the following figure 4).


Figure. 4

4. Operating handle component extends the loop and gently wraps the target tissue or polyps to be resection. Do not overexert the trap tissue or polyps. This can lead to injury, such as puncture wounds, bleeding, or injuries caused by mechanical removal. (as shown in the figure below)


Figure. 5

5. Confirm the output parameters of the high-frequency electrogenerator. If normal, the high-frequency electrogenerator can be started to remove the polyps or tissues.
6. After the resection, close the high-frequency electrogenerator and disconnect the electrode plug of the snare and the electrode socket of the high-frequency electrogenerator.
7. Put the loop into the outer tube and draw out the snare.

EMC CONDITIONAL

The Sterile Hot Snare meets emission standards of Class A and Group1.

Guidance and manufacturer's declaration - electromagnetic emissions –
for all ME EQUIPMENT and ME SYSTEM.

Table 1 – Emission limits per environment		
Phenomenon	Compliance	Electromagnetic environment
Conducted and Radiated RF	CISPR 11	Professional healthcare facility environment

EMISSIONS	Group 1, Class A	
Harmonic distortion	IEC 61000-3-2 Class A	
Voltage fluctuations and flicker	IEC 61000-3-3 Complies	

Guidance and manufacturer's declaration - electromagnetic immunity –
for all ME EQUIPMENT and ME SYSTEM.

Table 2 – Electromagnetic immunity		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM fields	IEC 61000-4-3	3V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz
Electrical fast transients/bursts	IEC 61000-4-4	±2kV, 100kHz (AC power port) ±1kV, 100kHz (Signal input/output parts port)
Surges	IEC 61000-4-5	±0.5kV, ±1kV (line to line); ±0.5kV, ±1kV, ±2kV(line to ground)
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz; 6V in ISM bands between 0.15MHz and 80MHz 80% AM at 1kHz
Voltage interruptions	IEC 61000-4-11	0% U _T ; 250/300 cycle
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz; 6V in ISM bands between 0.15MHz and 80MHz 80% AM at 1kHz
Voltage dips	IEC 61000-4-11	0% U _T ; 0.5cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°; 0% U _T ; 1cycle and 70% U _T ; 25/30 cycles Single phase: at 0°

Guidance and manufacturer's declaration – electromagnetic immunity for
ME EQUIPMENT and ME SYSTEM that are not LIFE-SUPPORTING.

Table 3 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment						
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation ^{b)} 18Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM ^{c)} ±5kHz deviation 1kHz sine	2	0,3	28

Table 3 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
710	704-787	LTE Band 13,17	Pulse modulation ^{b)} 217Hz	0,2	0,3	9
745						
780						
810	800-960	GSM800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18Hz	2	0,3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4,25;UMTS	Pulse modulation ^{b)} 217Hz	2	0,3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217Hz	2	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217Hz	0,2	0,3	9
5500						
5785						
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT of ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a) For some services, only the uplink frequencies are included.						
b) The carrier shall be modulated using a 50% duty cycle square wave signal.						
c) As an alternative to FM modulation, 50% pulse modulation at 18Hz may be used because while it does not represent actual modulation, it would be worst case.						

PRODUCT DISPOSAL

After use, this device may be a potential biohazard. Handle and dispose of in accordance with hospital, local and administrative, laws and regulations.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Do not re-use		Do not re-sterilize
	Date of Manufacture		Manufacturer
	Use-by date		Caution
	Catalogue number		Batch code
	This device is not made with natural rubber latex		Sterilized using ethylene oxide
	Diameter		Working length
	Compatible working channel		Authorized representative in the European Community/European Union
	Keep away from sunlight		Keep dry
	Do not use if package is damaged and consult instructions for use		Medical Device
	Sterile barrier system/sterile packaging		Open here
	Consult instructions for use or consult electronic instructions for use		Importer
	Type BF applied part		Unique device identifier
	Temperature limitation		Humidity limitation
	Atmospheric pressure limitation		Follow instructions for use
	Contents		

【Packaging】 Packed in Pouch

【Production Date】 See packaging

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 See packaging

WARRANTY

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the device is used by Buyer, the devices will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the devices and do not constitute any express warranties. Any technical advice with respect to the device and guarantee of specific properties of or in the devices shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for device failure or deficiency due to improper storage, alteration, or the consequences of uses for which the devices were not designed or that adversely affect the devices' integrity, reliability, or performance.

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Instruction page for electronic Instructions for Use (eIFU) website

To view / access IFUs, please go to <https://eifu.micro-tech.com.cn> by entering the product-specific IFU code HS01.

The hardware and software requirements to view / access IFUs are as follows:

- Both the computer or mobile devices (mobile phones, pad, notebooks and computer) can be used, with current operating system (computer: Apple's macOS system, Microsoft's Windows system; Mobile: Apple iOS system, Google Android system, HarmonyOS system).
- The browser must be the latest generation of Google Chrome, Firefox, Safari or Microsoft Edge. For a paper copy, provided in 7 days at no cost, please call +49(0)211 73 27 626-0.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp. GmbH (Europe)