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WARNINGS

1. The product is intended for single use only! DO NOT reuse, reprocess or resterilize the device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse or reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
2. Do not use this device for any purpose other than the stated indications.
3. Operation of this instrument is based on the assumption that open surgery is possible as an emergency measure if the clip cannot be detached from the instrument or if any other unexpected circumstance takes place.
4. It might be impossible to stop bleeding depending on the hemorrhage situation because the clip performance for hemostasis is limited. Prepare more than one hemostasis device.
5. Re-bleeding may occur on the clipping site, depending on the local condition. Check the patient for any re-bleeding after the operation as appropriate.
6. As the clips are made of stainless-steel, do not use them on a patient who is severely allergic to metals. Otherwise, allergic symptoms may occur.
7. This device does not contain latex.

DEVICE NAME

Sterile Repositionable Hemostasis Clipping Device

DEVICE DESCRIPTION**【CONTENTS】**

The Sterile Repositionable Hemostasis Clipping Device is a sterile, single-use endoscopic clipping device, intended to be used for hemostasis for mucosal/submucosal defects in digestive tract.

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The specifications of the product are represented by open width (L1), the maximum outside diameter (D) and working length, as indicated in table 1.

Applicable endoscope working channel: $\geq \varnothing 2.8\text{mm}$.

Table 1 Specification for Sterile repositionable hemostasis clipping device

Type			Open width (mm)	Maximal Outer Diameter of the Insertion Part D (mm)	Working Length (mm)
ROCC	High-end type Normal type Super Eco	C	8	2.6	1650
		D	11		1950
		F	16		2300
		G	20		2350
LOCK	High-end type Normal type Super Eco	C	8	2.6	1650
		D	11		1950
		F	16		2300
		G	20		2350

【STRUCTURE】

Sterile Repositionable Hemostasis Clipping Device mainly includes clip component, distal end of spring tube, proximal end of spring tube and handle component (Fig. 1).


Fig. 1 Schematic diagram of Sterile Repositionable Hemostasis Clipping Device

1. Clip Assembly 2. Distal end of spring tube 3. Proximal end of spring tube 4. Handle

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by trained health care professional only. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

【PERFORMANCE CHARACTERISTIC】

This device can be repositioned before deployment.

The clip assembly can follow the rotation when rotating the handle.

【MR SAFETY INFORMATION】

The vitro simulative test on single Sterile Repositionable Hemostasis Clipping Device has confirmed that, according to the ASTM F2503 standard, patients with one clip or multiple dispersive clips or even multiple centralized clips which not touch each other, could get safely scanned under magnetic resonance parameters as follows.

Static Magnetic Field

- Static magnetic field of 1.5-Tesla or 3-Tesla
- The maximum spatial gradient should be no more than 4000 G/cm.
- The maximum MRI system report: the averaged specific absorption rate in whole body level is 2 W/kg.

MRI-Related Heating:

After 15 minutes of continuous scan based on the parameters described as above, the increase of temperature caused by the clips is no more than 2°C.

Artifact Information:

In non-clinical testing, the image artifact caused by the device extends approximately 25mm from the clip when imaged with a gradient echo or spin echo pulse sequence in a 3 Tesla MRI system.

【RADIOPACITY INFORMATION】

The implanted part Clip Assembly of Sterile Repositionable Hemostasis Clipping Device in the human body is made of stainless steel and is visualized under X-ray, the visibility under X-ray can be clearly identified.

INTENDED USE /INTENDED PURPOSE

This device is intended to be used for hemostasis for mucosal/submucosal defects in digestive tract.

INDICATIONS

- (1) Hemostasis for
 - a. Mucosal / sub-mucosal defects < 3cm
 - b. Bleeding ulcers
 - c. Polyps < 1.5cm in diameter
- (2) As a supplementary method, closure of GI tract luminal perforations <20mm that can be treated conservatively

CONTRAINDICATIONS

1. Mucosal/submucosal defects greater than 3cm.
2. Polyps greater than 1.5cm in diameter.
3. The endoscopy cannot be tolerated or the endoscope cannot reach the target site(s).
4. Do not use this device when hemostasis cannot be verified visually with an endoscopic field of view.

ADVERSE EVENTS

The potential adverse effects associated with Sterile Repositionable Hemostasis Clipping Device may include:

Bleeding

Perforation

Pain

Infection

Septicemia

Inflammation of tissue

Mucosal damage

Complications which are not currently known or observed may be present.

PATIENT POPULATION

The product is only intended for adult and adolescent populations.

CLINICAL BENEFIT

The device can close perforation or close the wound to achieve hemostasis.

PRODUCT MATERIAL SUBSTANCE TABLE

Part Name	Overall Materials	Hazardous substance to which patients can be exposed		
		Bisphenol A	Phthalates	CMR&ED
Clip Assembly	SUS631,SUS316,SUS304	○	○	○
Delivery system (spring tube and handle)	SUS304, SUS316, SUS304+HDPE,ABS,ABS+TPE	○	○	○
Remarks: ○ = Not Detected or less than method limit or meet the limit requirement ● Phthalates and Bisphenol A meet the requirements of REACH EC /1907/2006, and the content of all substances shall not exceed 0.1%. ● CMR&ED substances meet the requirement of Medical Device Regulation 2017/745, and the content of CMR&ED substances in concentration is not above 0.1% weight by weight (w/w).				

OVERALL QUALITATIVE AND QUANTITATIVE INFORMATION

The overall qualitative and quantitative material information for the implanted clip assembly is as follows:

Implanted Part Name	Material	Content
Clip Assembly-LOCK	SUS631	30%-40%
	SUS304	1%-5%
	SUS316	60%-65%
Clip Assembly-ROCC	SUS631	25%-35%
	SUS304	2%-5%
	SUS316	65%-70%

CAUTIONS

1. This instruction manual contains essential information on using this instrument safely and effectively. Before use thoroughly review this manual and the manuals of all equipment which will be used during the procedure and use the instruments as instructed. Keep this and all related instruction manuals in a safe, accessible location. If you have any questions or comments about any information in this manual, please contact Micro-Tech.
2. Patients should be informed of the potential risks and complications, which may lead to injury, illness or death of the patient
3. Please read this instruction manual entirely before use.
4. Do not use this instrument when hemostasis cannot be verified visually within the endoscopic field of view after application.

1. Please check the shelf life before use. Do not use any device if the shelf life is expired.
2. The product is intended for single use only and one patient only. Do not reuse the device, nor share it by two or more patients.
3. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State.
4. This device is intended for use by trained health care professional only. As an accessory of digestive endoscopy, the product shall be used by professional familiar with the operation technique of digestive endoscopy.
5. Do not operate the spring tube and clip with excessive force in case of any damage to the device.
6. Always observe the endoscopic image during operation. Make sure that the device is present on the image and the operation is normal, If the clip falls off from the endoscope before release, take it out with foreign body forceps.
7. Please confirm that the visual field of endoscope is clear before use.
8. Dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
9. Limited studies indicate that clipping hard or severely fibrotic lesions to achieve hemostasis may be more difficult.
10. The number of clips required for hemostasis may vary depending upon the anatomical site histology, lesion type and patient condition and history. A sufficient quantity of clips should be prepared in consideration of all of these factors prior to the procedure.
11. The use of clips in the presence of bacterial contamination may potentiate or prolong infection.
12. Patients who use anticoagulant drugs before surgery should be used with caution.
13. A single copy of the instructions for use is provided in a case box (10 devices). If further copies are requested, please contact Micro-Tech we will provide free of charge.
14. If you need more clinical information about the product, you can view the Summary of Safety and Clinical Performance (SSCP). It will be uploaded to the European Database on Medical Devices (Eudamed) (link: <https://ec.europa.eu/tools/eudamed>), where it is linked to the Basic UDI-DI.

HOW SUPPLIED

Product packaging intact and has been sterilized (EO sterilization)

One-time use, for one patient to use only

STORAGE

1. Store the device in a cool, dark and dry place, and do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation. The sterilization is valid for 3 years. Prior to use, check the expiry date of sterilization on the label of the package.
2. Storage conditions: in a room with a temperature range of 10-25°C and a humidity range of 10%-80%.

COMPATIBILITY

The proposed device is intended to use in endoscopic surgery of digestive tract by surgeon. It is used in combination with digestive tract endoscope. The specifications of digestive tract endoscopes should comply with the following requirements:

- 1) Channel diameter $\geq 2.8\text{mm}$,

- 2) Working length < the working length of proposed device.

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
1. Check whether there is any damage on the packaging before use. Do not use the device if the packaging is damaged.
2. After the shelf life is confirmed, unpack and take the device out of the package.
3. Do not use the device if there are any function abnormalities or questions.



【INSTRUCTIONS FOR USE】

1. The device is compatible with an endoscope channel of 2.8mm or larger.
2. Remove protective sleeve. Carefully insert device into endoscope channel, ensuring that the clip is in the closed position (See Fig.2).


NOTE: Endoscope must remain as straight as possible when inserting the device.

3. Pass the clip through the endoscope channel towards targeted site. Open the clip pushing slider forward (See Fig.3).
4. Clip can be rotated clockwise or anti-clockwise by turning handle until correct position is achieved. Advance device until contact is made with targeted site (See Fig.4).


5. When satisfied with clip positioning, close the clip component against the tissue by pulling the slider back until tactile resistance is felt in the handle. Clip position may now be assessed prior to deployment. If clip is not in desired position, clip may be re-opened and repositioned (See Fig.5).

NOTE: Do not continue to pull back the slider further beyond the tactile resistance until you are ready to deploy the clip, otherwise you may not be able to re-open the clip. This is explained in the next part, point No.6.

6. To deploy the clip, continue pulling back on slider beyond the tactile resistance point. You will feel an audible snap as clip component detaches (See Fig.6).

NOTE1: If the clip did not immediately detach from the catheter, then apply gentle movement of the catheter or endoscope to unseat the clip.

7. Remove sheath from endoscope.

NOTE: Endoscope must remain as straight as possible when withdrawing the device.

PRODUCT DISPOSAL

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

POST-PROCEDURE

1. After operation, follow the doctor's advice not to eat food too early. Improper diet, forcible defecation, cough and other activities increasing abdominal pressure will result in forcible falling off of the clip and bleeding afterwards.
2. Patients are followed up regularly. Generally, some patients' clips will fall off spontaneously within 1 month and excreted out of the body with the digestive system, while a few clips are excreted out of the body within 3 months; In rare cases, the duration is more than three months, and if the clip does not fall off more than six months, the doctor can remove the clip with a grasping forcep. If the clip is removed with a grasping forcep, there may be a risk of bleeding. It is recommended that during the removal process, first try to gently shake the grasping forcep to loosen the adhesion site, and then slowly separate and remove it, minimizing bleeding caused by tissue damage. Continuously monitor the patient's vital signs after the operation and handle any abnormalities in a timely manner.

SYMBOL INDICATIONS

	Do not re-use		Do not re-sterilize
	Date of manufacture		Manufacturer
	Use-by date		Caution
	Catalogue number		Batch code
	MR conditional		Sterilized using ethylene oxide
	Authorized representative in the European Community/European Union		Keep away from sunlight
	Keep dry		Do not use if package is damaged
	Consult instructions for use or consult electronic instructions for use		Channel diameter
	Working length		Diameter
	This Product is not made with Natural Rubber Latex		Medical Device
	Sterile barrier system/ Sterile packaging		Importer
	Unique Device Identifier		Temperature limit: 10-25° C
	Humidity limitation: 10%-80%		Open here

【Packaging】 Flexible peel pouch

【Production Date】 See packing

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 Three years

WARRANTY

Micro-Tech Nanjing (MT) warrants that reasonable care and prudence has been exercised in the design and manufacture of this device. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning, and sterilization of this device as well as factors related to the patient, diagnosis, treatment, surgical procedures, and other matters beyond MT's control directly affect the device and the results obtained from its use. MT's obligation under this warranty is limited to the repair or replacement of this device and MT's shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this device. MT neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this device. MT assumes no liability with respect to devices reused, reprocessed or resterilized and makes no warranties, expressed or implied, including (but not limited to) merchantability or fitness for intended use, with respect to such device.

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Instruction page for electronic Instructions for Use (eIFU) website

To view / access IFUs, please go to <https://eifu.micro-tech.com.cn> by entering the product-specific IFU code ROCC03.

The hardware and software requirements to view / access IFUs are as follows:

- Both the computer or mobile devices (mobile phones, pad, notebooks and computer) can be used, with current operating system (computer: Apple's macOS system, Microsoft's Windows system; Mobile: Apple iOS system, Google Android system, HarmonyOS system).

- The browser must be the latest generation of Google Chrome, Firefox, Safari or Microsoft Edge.

For a paper copy, provided in 7 days at no cost, please call + 49(0)211 73 27 626-0.

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Micro-Tech (Nanjing) Co.,Ltd.
Micro-Tech Europe GmbH
Shanghai International Holding Corp. GmbH (Europe)