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**WARNING:**

---

1. **The product is intended for single use only!** DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization, or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness, or death. Re-use, re-sterilization, or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness, or death of the patient. Micro-Tech assumes no liability with respect to devices reused, re-sterilized, or reprocessed.
2. Do not use this device for any purpose other than its intended use.
3. Do not insert the device into the endoscope unless there is a clear endoscopic field of view. If the distal end of the insertion portion in the endoscopic field of view cannot be seen, do not use the device. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, bleeding, or mucous membrane damage. It may also damage the endoscope and/or device.
4. This device is not intended to be used in the presence of flammable liquids, in an oxygen enriched atmosphere, or in the presence of explosive gases. Any electrosurgical device constitutes a potential electrical hazard to the patient and/or operator.
5. Avoid High Frequency (HF) output settings where the maximum output voltage exceeds rated accessory voltage.
6. Use this device in an environment equipped to accommodate open surgery and have a hospitalization plan prepared in case a problem occurs that cannot be resolved endoscopically.
7. This device does not contain any user-serviceable parts. Do not disassemble, modify, or attempt to repair the device; patient or user injury and/or equipment damage may occur.
8. Do not force the cutting knife against tissue with excessive force while activating output. Otherwise, unintended resection, perforation, bleeding, and cracking of the device may occur. When resecting tissue, always confirm the direction of resection and use the device without excessive force.
9. If the instrument is used on a patient with an implanted pacemaker, serious harm to the patient may occur. This instrument may cause an implanted pacemaker to malfunction. Before proceeding, always confirm with a cardiologist or the manufacturer of the pacemaker that is safe.
10. When using the device in the vicinity of the heart, be sure to use it with the minimum necessary output. Spark discharge during operation may affect the heart.
11. When using an electrocardiograph or other physiological monitoring equipment simultaneously with the device on a patient, any monitoring electrodes should be placed as far away as possible from the electrodes used with the electrosurgical unit. Needle monitoring electrodes should not be used, as they may cause patient burns. Physiological monitoring equipment incorporating high-frequency current limiting devices is recommended.

**DEVICE NAME:**

---

Single Use Electrosurgical Knife

## DEVICE DESCRIPTION

### 【SPECIFICATIONS】

Table 1 Specification and parameter

Unit: mm

Specification	Cutting Knife Shape	Cutting Knife Length L1	Effective Working Length L	Injection capable or NOT	The maximum diameter of the insertion portion D
MK-I-1-165	I	1.5	1650	Yes	<2.7
MK-I-1-195	I	1.5	1950	Yes	<2.7
MK-I-1-235	I	1.5	2350	Yes	<2.7
MK-I-2-165	I	2	1650	Yes	<2.7
MK-I-2-195	I	2	1950	Yes	<2.7
MK-I-2-235	I	2	2350	Yes	<2.7
MK-I-4-165	I	4	1650	Yes	<2.7
MK-I-4-195	I	4	1950	Yes	<2.7
MK-I-4-235	I	4	2350	Yes	<2.7
MK-T-1-165	T	1.5	1650	Yes	<2.7
MK-T-1-195	T	1.5	1950	Yes	<2.7
MK-T-1-235	T	1.5	2350	Yes	<2.7
MK-T-2-165	T	2	1650	Yes	<2.7
MK-T-2-195	T	2	1950	Yes	<2.7
MK-T-2-235	T	2	2350	Yes	<2.7
MK-T-4-165	T	4	1650	Yes	<2.7
MK-T-4-195	T	4	1950	Yes	<2.7
MK-T-4-235	T	4	2350	Yes	<2.7
MK-O-1-165	O	1.5	1650	Yes	<2.7
MK-O-1-195	O	1.5	1950	Yes	<2.7
MK-O-1-235	O	1.5	2350	Yes	<2.7
MK-O-2-165	O	2	1650	Yes	<2.7
MK-O-2-195	O	2	1950	Yes	<2.7
MK-O-2-235	O	2	2350	Yes	<2.7
MK-O-4-165	O	4	1650	Yes	<2.7
MK-O-4-195	O	4	1950	Yes	<2.7
MK-O-4-235	O	4	2350	Yes	<2.7
MK-IT-1-165	IT	1.5	1650	Yes	<2.7
MK-IT-1-195	IT	1.5	1950	Yes	<2.7
MK-IT-1-235	IT	1.5	2350	Yes	<2.7

Specification	Cutting Knife Shape	Cutting Knife Length L1	Effective Working Length L	Injection capable or NOT	The maximum diameter of the insertion portion D
MK-IT-2-165	IT	2	1650	Yes	<2.7
MK-IT-2-195	IT	2	1950	Yes	<2.7
MK-IT-2-235	IT	2	2350	Yes	<2.7
MK-IT-4-165	IT	4	1650	Yes	<2.7
MK-IT-4-195	IT	4	1950	Yes	<2.7
MK-IT-4-235	IT	4	2350	Yes	<2.7
MK-I-1-165-N	I	1.5	1650	No	<2.7
MK-I-1-195-N	I	1.5	1950	No	<2.7
MK-I-1-235-N	I	1.5	2350	No	<2.7
MK-I-2-165-N	I	2	1650	No	<2.7
MK-I-2-195-N	I	2	1950	No	<2.7
MK-I-2-235-N	I	2	2350	No	<2.7
MK-I-4-165-N	I	4	1650	No	<2.7
MK-I-4-195-N	I	4	1950	No	<2.7
MK-I-4-235-N	I	4	2350	No	<2.7
MK-T-1-165-N	T	1.5	1650	No	<2.7
MK-T-1-195-N	T	1.5	1950	No	<2.7
MK-T-1-235-N	T	1.5	2350	No	<2.7
MK-T-2-165-N	T	2	1650	No	<2.7
MK-T-2-195-N	T	2	1950	No	<2.7
MK-T-2-235-N	T	2	2350	No	<2.7
MK-T-4-165-N	T	4	1650	No	<2.7
MK-T-4-195-N	T	4	1950	No	<2.7
MK-T-4-235-N	T	4	2350	No	<2.7
MK-O-1-165-N	O	1.5	1650	No	<2.7
MK-O-1-195-N	O	1.5	1950	No	<2.7
MK-O-1-235-N	O	1.5	2350	No	<2.7
MK-O-2-165-N	O	2	1650	No	<2.7
MK-O-2-195-N	O	2	1950	No	<2.7
MK-O-2-235-N	O	2	2350	No	<2.7
MK-O-4-165-N	O	4	1650	No	<2.7
MK-O-4-195-N	O	4	1950	No	<2.7
MK-O-4-235-N	O	4	2350	No	<2.7
MK-IT-1-165-N	IT	1.5	1650	No	<2.7
MK-IT-1-195-N	IT	1.5	1950	No	<2.7

Specification	Cutting Knife Shape	Cutting Knife Length L1	Effective Working Length L	Injection capable or NOT	The maximum diameter of the insertion portion D
MK-IT-1-235-N	IT	1.5	2350	No	<2.7
MK-IT-2-165-N	IT	2	1650	No	<2.7
MK-IT-2-195-N	IT	2	1950	No	<2.7
MK-IT-2-235-N	IT	2	2350	No	<2.7
MK-IT-4-165-N	IT	4	1650	No	<2.7
MK-IT-4-195-N	IT	4	1950	No	<2.7
MK-IT-4-235-N	IT	4	2350	No	<2.7

## 【STRUCTURE】

The Single Use Electrosurgical Knife includes a cutting knife assembly, an outer tube assembly, and a handle assembly (Fig.1 and Fig. 2).


Fig.1 Schematic diagram of Single Use Electrosurgical Knife (with injection function)

1. Cutting knife assembly 2. Outer tube assembly 3. Handle assembly


Fig.2 Schematic diagram of Single Use Electrosurgical Knife (without injection function)

1. Cutting knife assembly 2. Outer tube assembly 3. Handle assembly


Fig.3 Schematic diagram of cutting knife assembly

## 【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by healthcare professional who received appropriate training in digestive endoscopy or with her/his supervision. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

## INTENDED USE

---

These devices have been designed to be used with endoscopes and electrosurgical units for marking, dissection, elevation, irrigation and preparation of tissue layers in combination with monopolar cutting and coagulation within the digestive tract.

## INDICATIONS

---

The device is indicated to be used for the endoscopic submucosal dissection (ESD) and endoscopic mucosal resection (EMR) within the digestive tract.

## PATIENT TARGET GROUP

---

The device is intended to be use for adults.

## CONTRAINDICATIONS

---

Contraindications for these devices are the same as the endoscopic submucosal dissection (ESD) and endoscopic mucosal resection (EMR).

- 1) Severe coagulopathy
- 2) Severe cardiopulmonary disease
- 3) Those who cannot tolerate endoscopic examination
- 4) Any other factors that the doctor considers unsuitable for operation

## COMPLICATIONS

---

Infection, perforation, stricture, bleeding, mucosal damage, pain.

Complications which are not currently known or observed may be present.

## **CAUTIONS:**

---

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
3. Please read the instructions for use entirely before use.

## **HOW SUPPLIED**

---

The Single Use Electrosurgical Knife is supplied STERILE.

## **STORAGE**

---

The product should be stored in a clean, well-ventilated, non-corrosive gas environment. Do not store them in direct sunlight.

Do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

The product is valid for 3 years.

## **ENVIRONMENT**

---

Operation environment

Temperature limitation: 10 °C ~40 °C

Humidity limitation: 30%~85%

Atmospheric pressure limitation: 800hPa ~1060hPa

Transportation and storage environment

Temperature limitation: -40 °C ~60 °C

Humidity limitation: 30%~85%

Atmospheric pressure limitation: 800hPa ~1060hPa

## **COMPATIBILITY**

---

### **【Applicable Electrosurgical Unit】**

Compatible with Erbe VIO® 200 D Electrosurgical Units and its matched monopolar cable 20192-117.


Fig.4 Schematic diagram of compatibility size with the monopolar cable

### 【Rated High-Frequency Voltage】

CUT: maximum 1200Vp (2400Vp-p)

DO NOT use higher repeated peak voltage than the Cut- 1200Vp (2400Vp-p) .

Power: maximum 80 watts

COAG: maximum 1200Vp (2400Vp-p)

DO NOT use higher repeated peak voltage than the Coag- 1200Vp (2400Vp-p) .

Power: maximum 80 watts

### 【Applicable Endoscopes】

Endoscope that are legal listed in EUROPE are recommended such as Olympus, Fujifilm and Pentax.

Compatible Endoscopes Working Channel inner diameter :  $\geq \Phi 2.8\text{mm}$

Compatible Endoscopes Working Length :  $<$  the effective working length of the device

### 【Syringe】

If the device is used to inject saline, a syringe (with 6% luer connector)  $\leq 20\text{mL}$  (20cc) will be used to inject the saline.


Fig.5

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1. Choose the appropriate specification referring to the compatible working channel in package label.
2. Contents supplied STERILE.
3. Inspect the package for damage before use. Do not use if package is damaged. Open the package carefully after verifying the period of validity.
4. Before use, remove the distal protection tube and ensure that there are no loose or disconnected parts. Ensure that there are no sharp edges on the device before inserting into the patient.
5. Do not use the device if there are any signs of damage. Do not attempt to repair a nonfunctional or damaged device.
6. Hold the device and form a loop in the insertion portion of approximately 20cm in diameter. Operate the slider and confirm that the cutting knife extends and retracts smoothly. If the cutting knife does not operate smoothly and as intended, do not use the device and replace with a new device.
7. When injecting a solution, such as saline, into the injection lumen, attach a syringe ( $\leq 20$  ml/20 cc) to the Luer. Inject the solution to check the injection lumen and remove any air within the lumen.

### 【INSTRUCTIONS FOR USE】

1. Attach the patient plate.
2. Pull the slider to retract the cutting knife into the outer sheath. With the cutting knife retracted, carefully insert the device into the biopsy valve of the endoscope. Advance the device until the distal end of the insertion portion appears within the endoscopic field of view.

**Warning:**

- 1) Do not insert the device into the endoscope if the cutting knife is not completely retracted into the outer sheath. The distal end of the insertion portion may extend from the distal end of the endoscope abruptly. This could cause patient injury, such as perforations, bleeding, or mucous membrane damage. It may also damage the endoscope and/or device.
- 2) Do not abruptly angulate the endoscope's bending section while the distal end of the insertion portion of the device is extended from the distal end of the endoscope. This could cause patient injury, such as perforations, bleeding, or mucous membrane damage. It may also damage the endoscope and/or device.

**Caution:**

When inserting the device into the endoscope, hold it close to the biopsy valve and keep it as straight as possible relative to the biopsy valve. Otherwise, the device could be damaged.

3. Insert the plug of the knife into the electrosurgical unit's plug until it clicks.
4. Push the slider to extend the cutting knife.
5. Set the power switch of the electrosurgical unit to ON and set an appropriate power.

**Warning:**

- 1) Be sure to check the output power of the electrosurgical unit before use. If the unit is used without the proper output setting, perforation, bleeding, or mucous membrane damage may occur.
- 2) Do not set the output value of the electrosurgical unit too high or too low. Also, do not allow the activation time to be too long or too short. Set the high-frequency output mode of the electrosurgical unit optimally according to the conditions of the tissue to be cut. An excessive or insufficient output value may result in perforation, bleeding, mucous membrane damage, or thermal injuries to the non-target tissue.
- 3) Do not bundle the cord with cables from other medical equipment (electrocardiograph, endoscopic video system, electrosurgical unit, etc.) High-frequency signals and spark discharge noise during cauterization may cause malfunctions in other medical equipment that could have an adverse effect on the patient. This could also cause an abnormality in the output from the electrosurgical unit and could cause patient injury, such as perforation, bleeding, or mucous membrane damage.
- 4) When an irregularity is detected during the use of this device, do not continue to use the electrosurgical knife as perforation, bleeding or mucous membrane damage may occur.
- 5) Ensure that the electrosurgical unit is supplying energy to the device when making an incision. Incision without electricity may result in patient injury, such as perforation, bleeding, or mucous membrane damage.

**Precautions:**

- 1) To avoid burning healthy tissue, do not activate output if the cutting knife is in contact with non-target tissue.
- 2) Do not activate output when the distal end of the endoscope is too close to or in contact with body cavity tissue. This could burn the tissue and/or damage the endoscope.
- 3) Do not activate output if any of the patient's skin surfaces are touching each other (a bare arm and the side of the chest, for example). This could burn the patient.
- 4) Do not activate output when patient is in contact with metal parts of operating table or

other units. This could burn the patient, operator, or assistant.

- 5) Avoid using excessive force when removing tissue attached to the cutting knife. Applying excessive force, such as vigorously scraping the cutting knife with tweezers or abruptly and repeatedly extending and retracting the cutting knife, may cause the cutting knife to break or crack the distal end.
  - 6) Cutting and coagulation for long periods of time will increase the likelihood of charred tissue build up on the cutting knife. If injection capability is lost ensure the cutting knife is retracted and remove the device from the endoscope. Remove charred tissue from cutting knife.
6. Hold the cutting knife against the tissue and activate the high-frequency current to mark, dissect, and prepare the tissue layers.

**If utilizing an injection model, perform the following:**

If elevation is needed, insert the cutting knife into the desired position and inject the appropriate solution such as saline to elevate the mucous membrane, thereby helping to precisely cut.

If irrigation is needed, locate the device near the actual or the potential bleeding site and inject an appropriate solution such as saline to clear the endoscope field of view.

**NOTE:**

- 1) Reference current output levels in combination with ERBE electro-surgical unit.

Operative technique method	Applied organ	Mode	Effect	Output level
Marking	Esophagus, Stomach	SOFT COAG	Effect 5	20-50W
	Colon	SOFT COAG	Effect 5	20W
Incision/ Dissection	Esophagus, Stomach	ENDO CUT Q Cutting Duration 2 Cutting Interval 2	Effect 3	-
	Colon	ENDO CUT Q Cutting Duration 3 Cutting Interval 3	Effect 2	-
Coagulation	Esophagus, Stomach	FORCED COAG	Effect 2	40-60 W
	Colon	FORCED COAG	Effect 2	40 W

- 2) The recommended settings provided in the table are standard current output levels, which are used in the most common case to the best knowledge of Micro-Tech. When you operate the electro-surgical unit, always set an appropriate out level according to the following conditions:
  - The condition of tissue to be cut or coagulation
  - The type/configuration/rated high-frequency voltage of the device that you use
  - The contact area (length) between the cutting knife and the tissue
  - Operational conditions like use of injection solution, and so on
  - Your therapeutic strategy (whether you put a priority on the prevention of bleeding or on limiting thermal injury to surrounding tissue).

Using of dredging needle:

- If injection becomes difficult or the injection capability is lost, ensure the cutting knife is fully retracted and remove the product from the endoscope.

- Remove any charred tissue from the cutting knife and the distal end of the product surface. Then, with the cutting knife in its retracted state, insert it into the dredging needle hole. Move it back and forth several times to dredge the clogs.
- After dredging, inject saline with a syringe ( $\leq 20$  ml/20 cc) to flush the injection lumen until the injection function is restored.

Caution:

- 1) Prolonged cutting increases the likelihood of charred tissue buildup on the cutting knife tip. To minimize this, intermittently inject saline into the injection lumen during the activation process to reduce charred tissue accumulation at the knife tip.
- 2) Rough insertion or failure to properly align the distal end with the dredging needle before inserting the device may damage the needle stylet. Always ensure proper alignment and proceed slowly when inserting the distal end into the dredging needle.
- 3) Keep the cutting knife clean, as charred tissue buildup can reduce the device's effectiveness. Do not activate the device during cleaning to avoid injury to operating personnel.
- 4) If using this device to perform endoscopic submucosal dissection (ESD), inject an appropriate solution such as saline into the submucosal layer to elevate the mucous membrane prior to initiating cutting. If necessary, add injection of an appropriate solution such as saline during dissection. If dissection is performed while there is not enough space between the mucous membrane and muscle layer, perforations may occur.
- 5) Aspirate fluids, such as mucus, that adhere to the cutting knife, outer sheath, and body cavity tissues. Patient injury such as perforations, bleeding, mucous membrane damage, and thermal injury of tissue could result if output is activated when in contact with these adhering fluids. If current is discharged while the cutting knife is being separated from the mucosa under wet conditions, it may break the cutting knife or crack the distal end.
- 6) Always confirm that the knife and tissue to be cut are in the endoscopic field of view. Otherwise, perforations, bleeding, or mucous membrane damage may occur.
- 7) Only cauterize tissue after confirming that the cutting knife is in contact with the tissue.
- 8) Do not resect tissue too deeply. Deep resection of tissue may cause bleeding, perforation, pneumomediastinum, and/or aerodermection during or after procedure. When resecting tissue, confirm that there is no irregularity in the resecting area and monitor the patient's condition at all times.
- 9) Store unused active cutting knife in a location isolated from patient.

Warning:

- 1) Ensure that the position of the cutting knife is correct. If injection is in an inappropriate position or injection is performed too deep, patient injury may occur.
  - 2) Take appropriate measures if any fluids from the electro-surgical knife come in contact with the eye or skin of the operator, as patient fluids could pose an infection control risk and/or cause skin irritation.
7. Switch the electro-surgical unit to the OFF position. Pull the slider to retract the cutting knife into the outer sheath. With the cutting knife retracted, carefully withdraw the device from the endoscope.

Caution:

If the cutting knife is not retracted, it could damage the device.

**NOTE**

1. Do not coil the insertion portion with a diameter of less than 15 cm. This could damage the insertion portion.

2. When using simultaneously with accessories compatible with high-frequency currents, do not activate the output while the accessory is in contact with body cavity tissue or with this device. This may cause bleeding or thermal injury to the non-target tissue.
3. Do not forcibly push the slider when the cutting knife is extended from the outer sheath. It may damage the device.

## EMC CONDITIONAL

The Single Use Electrosurgical Knife meets emission standards of Class A and Group1.

Guidance and manufacturer's declaration - electromagnetic emissions –  
for all ME EQUIPMENT and ME SYSTEM.

<b>Table 1 – Emission limits per environment</b>		
Phenomenon	Compliance	Electromagnetic environment
Conducted and Radiated RF EMISSIONS	CISPR 11 Group 1, Class A	Professional healthcare facility environment
Harmonic distortion	IEC 61000-3-2 Class A	
Voltage fluctuations and flicker	IEC 61000-3-3 Complies	

Guidance and manufacturer's declaration - electromagnetic immunity –  
for all ME EQUIPMENT and ME SYSTEM.

<b>Table 2 – Electromagnetic immunity</b>		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM fields	IEC 61000-4-3	3V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz
Electrical fast transients/bursts	IEC 61000-4-4	±2kV, 100kHz (AC power port) ±1kV, 100kHz (Signal input/output parts port)
Surges	IEC 61000-4-5	±0.5kV, ±1kV (line to line); ±0.5kV, ±1kV, ±2kV(line to ground)
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz; 6V in ISM bands between 0.15MHz and 80MHz 80% AM at 1kHz
Voltage interruptions	IEC 61000-4-11	0% U <sub>T</sub> ; 250/300 cycle
Conducted disturbances induced by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz; 6V in ISM bands between 0.15MHz and 80MHz 80% AM at 1kHz

<b>Table 2 – Electromagnetic immunity</b>		
Phenomenon	Basic EMC standard	Immunity test levels
		Professional healthcare facility environment
Voltage dips	IEC 61000-4-11	0% $U_T$ ; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°; 0% $U_T$ ; 1 cycle and 70% $U_T$ ; 25/30 cycles Single phase: at 0°

Guidance and manufacturer's declaration – electromagnetic immunity for ME EQUIPMENT and ME SYSTEM that are not LIFE-SUPPORTING.

<b>Table 3 – Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment</b>						
Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380-390	TETRA 400	Pulse modulation <sup>b)</sup> 18Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM <sup>c)</sup> ±5kHz deviation 1kHz sine	2	0,3	28
710	704-787	LTE Band 13,17	Pulse modulation <sup>b)</sup> 217Hz	0,2	0,3	9
745						
780						
810	800-960	GSM800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation <sup>b)</sup> 18Hz	2	0,3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4,25;UMTS	Pulse modulation <sup>b)</sup> 217Hz	2	0,3	28
1845						
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation <sup>b)</sup> 217Hz	2	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation <sup>b)</sup> 217Hz	0,2	0,3	9
5500						
5785						

NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT of ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

**Table 3 – Test specifications for ENCLOSURE PORT IMMUNITY to  
RF wireless communications equipment**

Test frequency (MHz)	Band <sup>a)</sup> (MHz)	Service <sup>a)</sup>	Modulation <sup>b)</sup>	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
<p>a) For some services, only the uplink frequencies are included.</p> <p>b) The carrier shall be modulated using a 50% duty cycle square wave signal.</p> <p>c) As an alternative to FM modulation, 50% pulse modulation at 18Hz may be used because while it does not represent actual modulation, it would be worst case.</p>						

## **PRODUCT DISPOSAL**

---

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## **POST-PROCEDURE**

---

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

## SYMBOL INDICATIONS

	Do not re-use		Do not resterilize
	Date of Manufacture		Manufacturer
	Use-by date		This device is not made with natural rubber latex
	Catalogue Number		Batch code
	Sterilization using ethylene oxide		Authorized representative in the European Community/ European Union
	CE sign with the number of the notified body		Keep away from sunlight
	Keep dry		Do not use if package is damaged and consult instruction for use
	Type BF applied part		Compatible working channel
	Consult instructions for use or consult electronic instructions for use		Working length
	Caution		Contents
	Temperature limitation		Humidity limitation
	Atmospheric pressure limitation		Unique device identifier
	Sterile barrier system/ sterile packaging		Medical Device
	Open here		Importer

**【Packaging】** Flexible peel pouch

**【Production Date】** See packaging

**【Sterilization】** Sterilized by EO (ethylene oxide) gas

**【Period of Validity】** 3 years

**WARRANTY:**

---

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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## Instruction page for electronic Instructions for Use (eIFU) website

To view / access IFUs, please go to <https://eifu.micro-tech.com.cn> by entering the product-specific IFU code MK02.

The hardware and software requirements to view / access IFUs are as follows:

- Both the computer or mobile devices (mobile phones, pad, notebooks and computer) can be used, with current operating system (computer: Apple's macOS system, Microsoft's Windows system; Mobile: Apple iOS system, Google Android system, HarmonyOS system).
  - The browser must be the latest generation of Google Chrome, Firefox, Safari or Microsoft Edge.
- For a paper copy, provided in 7 days at no cost, please call + 49(0)211 73 27 626-0.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp. GmbH (Europe)