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## WARNING

---

1. The product is intended for single use only! DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed.
2. Do not use this device for any purpose other than its intended use.
3. Do not insert the device into the endoscope unless you have a clear endoscopic field of view. If you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhages or mucous membrane damage. It may also damage the endoscope and/or device.
4. The operation of HF Generator may adversely influence the operation of other electronic equipment. For PATIENTS with cardiac pacemakers or other active implants, a possible HAZARD exists because interference with the action of the active implant may occur. or the active implant may be damaged. In case of doubt, qualified advice should be obtained. Cardiologist should be consulted prior to using device in these patients.

## DEVICE NAME

---


Sphincterotome / short-wire compatible

## DEVICE DESCRIPTION

---

### 【MODEL NUMBER】


RDSP-□□ □□-□ □ □

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Cutting wire			Outer tube			Preloaded Guidewire Length & Diameter
Length (±2 mm)	Cutting wire type	Insulation tube	Diameter (±0.2 mm)	Working length (±100mm)	Tip length (+2/0 mm)	
20, 25, 30	Single, Weave	With insulation tube, Without insulation tube	2.5	2000	0, 2, 5, 10	No, ≤0.035inch&2600mm, ≤0.035inch&4500mm

**【STRUCTURE】**

DETAIL I

SCALE 15:1

Figure 1 Product Drawing of Sphincterotome / short-wire compatible

No	Component	No	Component	No	Component
1	Cutting wire	6	Injection luer port	11	X-ray developing ring
2	Insulation tube	7	Limited block	12	Connecting pipe
3	Outer tube	8	Electric conduction seat plug	13	Guidewire
4	Heat shrink tube	9	Finger loop	14	Mandrel
5	Guidewire port	10	Core bar	15	/

**INTENDED USE/INDICATIONS FOR USE**

The device is intended to be used with endoscope and guidewire for selective cannulation of the biliary ducts and monopolar cutting in sphincterotomy of the Papilla of Vater and/or the Sphincter of Oddi using high-frequency current. The device can also be used to inject contrast medium.

## CONTRAINDICATIONS

---

Contraindications for this device are those for endoscopic retrograde cholangiopancreatography (ERCP) and endoscopic sphincterotomy (ES), include:

1. Patients who have obstruction of upper digestive tract, the device may not reach the descending duodenum;
2. Inability to pass the endoscope or tolerate the endoscopy;
3. Severe cardiopulmonary disease;
4. Severe coagulopathy;
5. Acute pancreatitis and severe biliary tract infection.
6. Severe allergic reaction to contrast medium;
7. Any other conditions the physician judges unsuitable for use.

## ADVERSE EVENTS

---

- ◆ Pancreatitis
- ◆ Perforation
- ◆ Hemorrhage
- ◆ Cholangitis
- ◆ Sepsis
- ◆ Hyperamylasemia
- ◆ Burns
- ◆ Stimulation
- ◆ Arrhythmias
- ◆ Complications which are not currently known or observed may be present

## CAUTION

---

1. Please read the instructions for use entirely before use.
2. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
3. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
4. Do not modify this device in any way, form or matter.
5. This device is not intended to be used in the presence of Flammable liquids, in an Oxygen Enriched Atmosphere or in the presence of explosive gases. Any electrosurgical device constitutes a potential electrical hazard to the patient and/or operator.
6. Avoid High Frequency output settings where maximum output voltage exceed rated accessory voltage.
7. The device is intended for adult and adolescent populations.
8. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

## HOW SUPPLIED

---

The device is supplied sterilized.

## COMPATIBILITY

---

### **【Applicable Electrosurgical Unit】**

High-frequency electrosurgical units that are legal marketed in EUROPE are recommended, such as ERBE VIO® 200 D.

### **Rated Accessory Voltage of High-frequency electrosurgical unit:**

CUT:1200Vp (2400Vp-p)

DO NOT use repeated peak voltage higher than the Cut- 1200Vp (2400Vp-p)

Power: 30W~70 W

### **【Matched monopolar cables model】**

ERBE 20192-117, CONMED 474-L.

Please check whether the neutral electrode used is compatible with the VIO High-frequency electrosurgical units.

### **【Applicable Endoscopes】**

Endoscope which is legal listing in European Union is recommended such as Olympus.

Compatible Working Channel:  $\geq \Phi 3.2\text{mm}$

Working Length: < the working length of The device

### **【Compatible Guidewire】**

The diameter of guidewire should be no more than 0.035inch(0.89mm), also the indication of the guidewire should be for biliary tract system.

Remark: For the models with Preloaded guidewire, the guidewire is 0.035inch, no additional guidewire will be used during the clinical use.

For the models without preloaded guidewire, a 0.035inch guidewire with the working length no less than 2600mm will be used during the clinical use.

### **【Guidewire Locking Device】**

Compatible with the endoscope

### **【Syringe】**

If the device is used to inject contrast medium, a syringe  $\leq 20\text{mL}$  (20cc) will be used to inject the contrast medium

## PRECAUTIONS

---

The device shall only be used by physicians or under the supervision of physicians trained in ERCP or ES. A thorough understanding of the technical principles, clinical applications and risks associated with ERCP / ES is necessary before using this device.

Safe and effective electrosurgery is dependent not only on equipment design, but also, to a large extent, on factors under the control of the operator. It is important that the following be read, understood, and followed in order to enhance safety and effectiveness:

1. Monitor sphincterotomy under fluoroscopy.
2. Avoid bowing of the device while it is in a coiled position.
3. Any use of this device, other than those listed in this instruction, is not recommended.
4. Recommend HF (high frequency) Generator model ERBE electrosurgical.
5. Inspect the device before use. If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify MT for return authorization.
6. The Rated Accessory Voltage for this device is 1200Vp(2400Vp-p).
7. This device is not intended to be used in the presence of Flammable liquids, in an Oxygen Enriched Atmosphere or in the presence of explosive gases. Any electrosurgical device constitutes a potential electrical hazard to the patient and/or operator.
8. This device complies with the EMC standard (IEC 60601-1-2:2014) for medical electrical equipment; When connecting to an instrument that complies with the EMC standard for medical electrical equipment, the whole system complies with the EMC standard for medical electrical equipment.
9. The output power selected should be as low as possible for the intended purpose. Certain devices or accessories may present an unacceptable RISK at low power settings. For example, with argon beam COAGULATION, the risk of gas embolism rises if there is insufficient HF power to produce a rapid, impermeable eschar on the target tissue.
10. Skin-to-skin contact should be avoided (for example between the patient's arms and body) by way of dry cloth or gauze. Monitoring electrodes should be placed as far from the surgical area as possible. Needle monitoring electrodes are not recommended. Avoid incidental contact between Active Cords and the patient's body, or any other electrodes. Flammable agents used for cleaning or disinfecting, or as solvents of adhesives, should be allowed to evaporate before the procedure.
11. Consult the neutral electrode manufacturer about the proper grounding of the patient. It is recommended that a monitoring neutral electrode be used, if a contact quality monitor is available, or built into the generator. The entire area of the neutral electrode should be attached reliably to the patient's body, and as close to the operating field as possible. The patient should not come into contact with metal parts or objects that may be grounded to earth. The use of antistatic sheeting is recommended for this purpose. Ensure that a good return path to the electrosurgical unit is maintained.
12. When HF Generator and physiological monitoring equipment are used simultaneously on the same PATIENT, any monitoring electrodes should be placed as far as possible from the surgical electrodes. Needle monitoring electrodes are not recommended.  
In all cases, monitoring systems incorporating HF current limiting devices are recommended.
13. The PATIENT leads should be positioned in such a way that contact with the PATIENT or other leads is avoided. Apparent low output or failure of the HF Generator to function correctly at the normal operating settings may indicate faulty application of the NEUTRAL ELECTRODE or poor contact in its connections. In this case, the application of the NEUTRAL ELECTRODE and its connections should be checked before

selecting a higher output power.

## PREPARATION

---

1. Inspect the package before use for any damage. Do not use if package is damaged. If this device shows any signs of damage, do not use. Do not attempt to repair a nonfunctional or damaged device.
2. Verify the expiration date. Do not use if expired.
3. Remove the device from the package.

Note: Ensure the pre-curved mandrel portion stays with the guidewire lumen of the device.

4. Remove the metal mandrel out of the device distal tip (Be careful to keep pre-curved shape of the tip)
5. Inspect the device for proper deflection and for any damage, such as kinks.  
Precaution: Kinks in the catheter will hinder injection capability. Do not use the device if any defect is found during inspection. Please notify manufacture and return for replacement.  
Precaution: Distal tip deflection greater than 90 degrees is not recommended.  
WARNING: Distal tip deflection greater than 90 degrees may cause the cutting wire to break.  
Precaution: Do not pre-active the device. Pre-activating the cutting wire prior to use may cause premature cutting wire fatigue and may compromise the cutting wire's integrity.

6. For the models without preloaded guidewire, a guidewire with diameter  $\leq 0.035$ " should be ready for use and be inserted it into the guidewire lumen of the device.
7. Ensure that the patient is properly grounded per the electrosurgical unit manufacturer's instruction. Recommend HF (high frequency) Generator model ERBE electrosurgical.
8. Connect a suitable electrosurgical active cord as required.
9. If contrast media will be injected through the injection luer port, attach a syringe  $\leq 20$  mL (20cc) to the proximal luer hub and inject some saline to check the injection luer port. Contrast medium which is legal in EUROPE is recommended such as iohexol.
10. The device is now ready to use.

## DIRECTIONS FOR USE

---

**Precaution:** It is recommended that the electrosurgical unit be turned off during insertion of the device. Advance the device using short deliberate 20mm-30mm movements to prevent kinking in the catheter.

**WARNING:** It is not allowed to open the foot switch during insertion of the device.

1. Ensure the Guidewire Locking Device is firmly seated onto the endoscope handle against the side of the working channel port; Insert the device into the endoscope with a working channel of 3.2mm or higher.
2. Direct the cutting wire position to the desired orientation by rotating the handle either left or right.
3. Once the device is at the desired position within the biliary system, use contrast medium to perform cholangiogram with fluoroscopy and/or x-ray or the guidewire to determine proper position within the biliary system.  
**WARNING:** Verify that the cutting wire has exited the endoscope by visualizing it on the endoscope monitor. Failure to do so may result in contact between the cutting wire and endoscope while electrical current is applied. This may cause a short circuit, which may result in patient injury, broken cutting wire, and/or damage to the endoscope.
4. Prior to performing device, make sure the device cutting wire is in proper position when viewed

endoscopically.

**Precaution:** If the guidewire is to be removed during sphincterotomy, turn off the power on the electro-surgical unit before removal.

- Using the appropriate power setting advised by the electro-surgical unit manufacturer (30W~70W), activate the cutting wire to perform sphincterotomy.

WARNING: Too low, maybe the cutting is not ideal, Too high, breakage risk may occur.

- Upon completion of the sphincterotomy, turn off the electro-surgical unit, disconnect the active cord, and then remove the device.
- To remove the device, separate the guidewire from the open channel on the wall of the device's catheter. Lock the guidewire on the Guidewire Locking Device, withdraw the device until the end of open channel emerges in the device entrance of Guidewire Locking Device, unlocking Guidewire and remove the device from endoscope completely
- If the device in use appeared wire broken, the head severely deformed or other damage, please do not use and replace the device. Discard the device according to local laws and/or regulations.
- If no other device exchange is performed, unlock the guidewire from the guidewire locking device and remove from the working channel of the endoscope.

**EMC CONDITIONAL**

Guidance and manufacturer's declaration - electromagnetic emissions - for all ME EQUIPMENT and ME SYSTEM.

**【WARNING】**

- The Operator should not use the system and should inform the customer service, if the monopolar cutting is lost or degraded due to EM disturbances.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 300mm to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Table 1 Electromagnetic emission level

<b>Electromagnetic emission</b>	
Electromagnetic requirements of this RF generator are given below and it is the responsibility of end user to meet these requirements.	
<b>Emission test</b>	<b>Compliance</b>
CISPR 11 Conducted emission	Group 1, Class A*
CISPR 11 Radiated emission	
IEC61000-3-2 Harmonic emission	N/A
IEC61000-3-3 Voltage fluctuation / flickering emission	N/A
* The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals	

(CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

Table 2 ENCLOSURE PORT

Phenomenon	Basic EMC standard	Immunity compliant levels
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF EM fields	IEC 61000-4-3	3 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See Table 4
RATED power frequency magnetic fields	IEC 61000-4-8	30 A/m 50 Hz or 60 Hz

Table 3 PATIENT coupling PORT

Phenomenon	Basic EMC standard	Immunity compliant levels
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz

Table 4 Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment

Test Frequency (MHz)	Band (MHz)	Service	Modulation	IMMUNITY TEST LEVEL (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 – 470	GMRS 460, FRS 460	FM ±5 kHz deviation 1 kHz sine	28
710	704 – 787	LTE Band 13,17	Pulse modulation 217 Hz	9
745				
780				

810	800 – 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3,4, 25; UMTS	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 – 5800	WLAN 802.11a/n	Pulse modulation 217 Hz	9
5500				
5785				

## STORAGE

Store in a clean and dry environment. Do not store in direct sunlight.

Transportation and storage environment	Ambient temperature	-35 °C~60 °C (-31 °F~140 °F)
	Relative humidity	10%~85%
	Atmospheric pressure	700hPa~1060hPa (0.7kgf/cm <sup>2</sup> ~1.1kgf/cm <sup>2</sup> ) (10.2psia ~15.4psia)
Operating environment	Ambient temperature	10 °C~40 °C (50 °F~104 °F)
	Relative humidity	30%~85%
	Atmospheric pressure	700hPa ~1060hPa (0.7kgf/cm <sup>2</sup> ~1.1kgf/cm <sup>2</sup> ) (10.2psia~15.4psia)

## PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

## SYMBOL INDICATION

	Manufacturer		Authorized representative in the European Community/ European Union
	Do not re-use		Do not re-sterilize
	Keep away from sunlight		Keep dry
	Use-by date		Date of manufacture
	Catalogue number		Batch code
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged and consult instructions for use
	Sterilized using ethylene oxide		Caution
	Temperature limit		Working Length
	Type BF applied part		Humidity limitation
	Sterile barrier system/ sterile packaging		Open here
	Medical Device		Importer
	Unique device identifier		Refer to instruction manual/booklet

**【Packaging】** Packed in flexible peel pouch

**【Production Date】** See packaging

**【Sterilization】** Sterilized by EO (ethylene oxide) gas

**【Period of Validity】** 3 Years

## **WARRANTY**

---

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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**Micro-Tech (Nanjing) Co., Ltd.**
**MICRO-TECH Europe GmbH**
**Shanghai International Holding Corp. GmbH (Europe)**