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WARNINGS

1. In accordance with the laws and regulations, the product is intended for use by trained professional doctor only.
2. The product is intended for single use only! DO NOT reuse, reprocess or resterilize the device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse or reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
3. Do not use this device for any purpose other than the stated indications
4. Operation of this instrument is based on the assumption that open surgery is possible as an emergency measure if the device cannot be detached from the instrument or if any other unexpected circumstance takes place.
5. Confirm that the endoscopy view is clear before use. Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. Damage to the endoscope and/or the instrument may also occur.
6. Carefully examine the unit to verify that neither the contents nor the package has been damaged in shipment. If the device shows any signs of damage DO NOT use this product and contact your local Micro-Tech representative.

DEVICE DESCRIPTION

[COMMON NAME]

Sterile Elastic Traction Device

DEVICE DESCRIPTION

[MODEL NUMBER]

ETD01006/ ETD01007/ ETD01008/ ETD01009/ ETD0101A/ ETD0101B


[SPECIFICATION]

Unit : mm

REF	Open width L1	Working Channel	Working Length L	Whether the Clip with Teeth	The center distance between the assembly ring and the last working ring
ETD01006	11±1	≥3.2	1650±50	With teeth	12.3±0.5
ETD01007	11±1	≥3.2	1950±50	With teeth	12.3±0.5
ETD01008	11±1	≥3.2	2350±50	With teeth	12.3±0.5
ETD01009	11±1	≥3.2	1650±50	Without teeth	12.3±0.5
ETD0101A	11±1	≥3.2	1950±50	Without teeth	12.3±0.5
ETD0101B	11±1	≥3.2	2350±50	Without teeth	12.3±0.5

[STRUCTURE]

The Product contains a traction clip and an elastic ring. Refer to Fig.1.(Fig. 2 is an enlarged view of the elastic ring)

L-Working Length D-Maximum Outer Diameter of Insertion Part L1 -Open Width

Fig.1 Schematic diagram of Sterile Elastic Traction Device


Fig.2 Elastic Ring

[USER INFORMATION/TRAINING/QUALIFICATIONS]

This device is intended for use by trained health care professional only. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

INTENDED USE

The Sterile Elastic Traction Device is indicated for use in endoscopic treatment to traction the tissue and expose the surgical view.

INDICATIONS

The Sterile Elastic Traction Device is indicated for use in flexible Endoscopy to provide retraction to assist in tissue resection, exposure, and removal of tissue within the gastrointestinal tract.

PATIENT TARGET GROUP

This device is intended to be used for adults.

CONTRAINDICATIONS

1. The patient with poor general condition cannot tolerate endoscopy.
2. The patient has narrow upper digestive tract where endoscope cannot pass through.
3. The patient has serious coagulation disorders and hemorrhagic diseases.
4. The patient is allergic to the device and the drugs used in the operation.
5. The patient who is not suitable to use the product according to the diagnosis.

ADVERSE EVENTS

The potential adverse effects associated with Sterile Elastic Traction Device may include:

Bleeding

Perforation

Pain

Infection

Septicemia

Inflammation of tissue

Mucosal damage

CAUTION

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
3. Please read the instructions for use entirely before use.

NOTES

1. The product is intended for use under the direct supervision of a suitably trained physician only. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Please check the shelf life before use. Do not use any device if the shelf life is expired.
3. The product is intended for single use only and one patient only. Do not reuse the device, nor share it by two or more patients.
4. Dispose of product and packaging in accordance with hospital, administrative and / or local government policy.
5. Do not operate the spring tube and clip with excessive force in case of any damage to the device.
6. Always observe the endoscopic image during operation. If the clip deploys prematurely, remove it with foreign body retrieval forceps.
7. For Usage in a side viewing endoscope: before actuating the elevator to the full position, assure that the clip body has projected adequately past the distal end of the scope, to prevent inadvertent detachment of the clip or damage to the device. Fully lower the elevator prior to removing the device following clip

deployment.

8. A notice to the user and / or patient that any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and / or patient is established.

HOW SUPPLIED

The product is supplied in sterile and intended for single use only.

STORAGE

The product should be stored in a cool, dry, clean, well-ventilated, non-corrosive gas environment.

COMPATIBILITY

The device is intended to use in endoscopic surgery of digestive tract and in conjunction with another Endoscopic Clip which should be able to pass through the endoscope channel smoothly and reach the target position. The specifications of digestive tract endoscopes should comply with the following requirements:

- 1) Channel diameter $\geq 3.2\text{mm}$.
- 2) Working length $<$ the working length of proposed device.

DIRECTIONS FOR USE

[PREPARATIVE CHECK AND PREPARATION]

Read the product instruction book and label, confirm the specification and the expiry date of sterilization. If the package is damaged or opened unintentionally before use, replace it with a new device.

[INSTRUCTIONS FOR USE]

1. Open the pouch and remove slider blister from the distal portion of the clip to expose the clip with the pre-loaded elastic ring. Carefully insert the device into the endoscope channel, ensuring that the clip is in the closed position (Fig.3).


Fig.3

Note: Wet the channel and distal portion of device with normal saline before inserting the device into the channel.

2. Advance the clip in small incremental movements towards the target site. Once in the instrument channel, there is no need to apply closure pressure on the handle.

Note: Don't retract the device in the inserting process.

Note: Applying excessive closure pressure to the handle during insertion may result in detachment of the clip.

Note: Endoscope should remain as straight as possible when inserting the device.

Note: When introducing the device in an endoscope in a tortuous position, straightening the endoscope may improve passage and exposure of the clip. With the clip in endoscopic view, carefully reposition the endoscope for treatment.

3. Endoscopically confirm that the device has exited the endoscope with the preloaded elastic ring still attached. If for some reason, the device exits the endoscope without the elastic ring attached, remove the device and retrieve the loop is appropriate. Select another device.

4. The clip can be rotated clockwise or counter-clockwise by slowly turning the handle component until the desired position is achieved. During rotation, the handle component and finger ring should be allowed to rotate.

5. Position distal tip of the device toward the targeted site.

6. Open the clip, with preloaded elastic ring attached and advance device into contact with the targeted site to properly anchor the clip to the tissue.

7. When satisfied with clip positioning, close the clip onto the tissue by pulling the finger rings back until tactile resistance is felt in the handle. The clip position may now be assessed prior to deployment (Fig.4). If clip is not in desired position, gently reopen the finger rings and repositioned the device.


Fig.4

Note: Do not continue to pull back the finger ring beyond the tactile resistance until you are ready to deploy the clip otherwise you may not be able to re-open the clip. If you hear or feel a click, the clip cannot be re-opened. To deploy the clip, continue pulling back the finger rings beyond the tactile resistance point. You will hear an audible snap when the clip component detaches.

Note: If separation of clip is not immediate, gently move catheter back and forth or use other endoscopic maneuvers necessary to unseat the catheter from clip.

8. At this point, the primary clip should be anchored to the desired. The preloaded, looped elastic ring is attached to one of the arms/jaws/projections of the anchored clip.

9. Open the pouch of another Endoscopic Clip. With clip closed, advance device in small increments into the accessory channel of the endoscope. The endoscopic is used to engage the elastic ring of the previously anchored clip.

Note: Please refer to the instructions for the operation method of Endoscopic clip.

10. Once the arm/jaw/projection of the Endoscopic clip has engaged with/traversed through the distal loop of the elastic ring it may be properly positioned to be engage into tissue to provide the desired traction/counter-traction/triangulation. Once the endoscopic clip is engaged with the targeted tissue and connected to the elastic ring of the primary clip, it is deployed from the catheter assembly as per step 7 above (Fig.5 and Fig.6). As the targeted tissue is cut and freed from its position the tension from the loop connecting the two clips, raises the flaps of tissue to expose the underlying tissue and proper cutting angles to assist in the resection process.


Fig 5


Fig.6

There are two traction schemes for reference:

Scheme 1: the elastic ring was used for "pathological auto traction"(Fig.7).

Scheme 2: the elastic ring was used for "traction from lesion site to normal site"(Fig.8).


Fig 7


Fig.8

Note: Do not over extend the loop. Over extending the loop may cause the loop to break and result in tissue damage or injury.

11. After the operation, the pathological tissue of scheme 1 can be removed according to the conventional ESD operation method. The pathological tissue of scheme 2 can be removed from the wall of digestive tract by foreign body forceps or snares (if bleeding occurs during operation, electrocoagulation or hemostasis clip), and then the pathological tissue and clips can be taken out together.

Note: Remove the clips carefully to avoid damaging the specimens before pretreatment.

PRODUCT DISPOSAL

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local laws and regulations.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Do not re-use		Do not resterilize
	Date of Manufacture		Manufacturer
	Use-by date		Not made with natural rubber latex.
	Contents		Batch code
	Catalogue Number		Sterilization using ethylene oxide
	Medical Device		Keep away from sunlight
	Authorized representative in the European Community		Do not use if package is damaged
	Keep dry		Compatible working channel
	Consult instructions for use or consult electronic instructions for use		Sterile barrier system
	Diameter		Caution
	Importer		Unique Device Identifier

[Packaging] Flexible peel pouch

[Production Date] See packaging

[Sterilization] Sterilized by EO (ethylene oxide) gas

[Period of Validity] 3 years

WARRANTY

Micro-Tech Nanjing (MT) warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness. Handling, storage, cleaning, and sterilization of this instrument as well as factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond MT's control directly affect the instrument and the results obtained from its use. MT's obligation under this warranty is limited to the repair or replacement of this instrument and MT's shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. MT neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. MT assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp.GmbH(Europe)