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WARNINGS

1. For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in inpatient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, reprocessed or resterilized.
2. This device shall not be used for any purpose other than its intended use.
3. Make sure the endoscope is clear before the instrument is inserted into the endoscope. Do not use the instrument if you cannot see the distal part of the instrument inserted in the endoscope. Unclear vision when the instrument is inserted into the endoscope may result in injury to the patient, such as perforation, hemorrhage or mucosal damage, as well as damage to the endoscope and/or instrument.
4. This device is intended for adult.
5. Patients should be informed of the potential risks and complications, which may lead to injury or illness of the patient.
6. The operation should be performed under good surveillance and necessary rescue equipment should be prepared in advance to ensure the safety and success of the operation.

DEVICE NAME

Extraction Basket

DEVICE DESCRIPTION**【CATALOGUE NUMBER】**

Table 1 Catalogue number

CEB00000	CEB00001	CEB00002	CEB00003	CEB00004	CEB00005	CEB00006
CEB01010	CEB01011	CEB01012	CEB01013	CEB01014	CEB01015	CEB01016
CEB02020	CEB02021	CEB02022	CEB02023	CEB02024	CEB02025	CEB02026
CEB03030	CEB03031	CEB03032	CEB03033	CEB03034	CEB03035	CEB03036
CEB04040	CEB04041	CEB04042	CEB04043	CEB04044	CEB04045	CEB04046

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The specifications of the product are represented by the Basket Diameter and Working Length, as indicated in Table 2.

Table 2 Specification and parameter

REF	Basket Shape	Basket Diameter (mm)	Basket Length (mm)	Compatible Endoscopes Working Channel Diameter (mm)	Working Length (mm)
CEB00000	Four wire spiral	10	30	≥1.1	2900
CEB01010		15	40		
CEB02020		20	50		
CEB03030		25	60		
CEB04040		30	70		
CEB00001		10	30		2000
CEB01011		15	40		
CEB02021		20	50		
CEB03031		25	60		

REF	Basket Shape	Basket Diameter (mm)	Basket Length (mm)	Compatible Endoscopes Working Channel Diameter (mm)	Working Length (mm)
CEB04041		30	70		1600
CEB00002		10	30		
CEB01012		15	40		
CEB02022		20	50		
CEB03032		25	60		
CEB04042		30	70		
CEB00003		10	30		1200
CEB01013		15	40		
CEB02023		20	50		
CEB03033		25	60		
CEB04043		30	70		
CEB00004		10	30		1000
CEB01014		15	40		
CEB02024		20	50		
CEB03034		25	60		
CEB04044		30	70		
CEB00005		10	30		800
CEB01015		15	40		
CEB02025		20	50		
CEB03035		25	60		
CEB04045		30	70		

REF	Basket Shape	Basket Diameter (mm)	Basket Length (mm)	Compatible Endoscopes Working Channel Diameter (mm)	Working Length (mm)
CEB00006		10	30		600
CEB01016		15	40		
CEB02026		20	50		
CEB03036		25	60		
CEB04046		30	70		

【STRUCTURE】

The Extraction Basket is composed with Handle assembly, Steel wire, Outer tube assembly and Basket assembly (See Figure 1 and Figure 2).


Figure 1 Schematic diagram of Extraction Basket

1. Handle assembly; 2. Steel wire; 3. Outer tube assembly; 4. Basket assembly


Figure 2 Schematic diagram of basket shape: Four wires spiral

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by physicians or medical personnel who have received appropriate training, in accordance with the applicable laws. As a consumable of digestive endoscope, the product shall be used by professionals familiar with the operation technique of digestive endoscope.

INTENDED USE/INDICATIONS

This device is used for the endoscopic removal of stones in the biliary system and foreign bodies.

PATIENT TARGET GROUP

This device is used for adult.

CONTRAINDICATIONS

1. Patients who have obstruction of upper digestive tract, the device may not reach the descending duodenum;
2. Inability to pass the endoscope or tolerate the endoscopy;
3. Severe cardiopulmonary disease;
4. Severe coagulopathy;
5. Acute pancreatitis and severe biliary tract infection;
6. Any other conditions the physician judges unsuitable for use.

COMPLICATIONS

1. Pancreatitis
2. Perforation
3. Hemorrhage
4. Cholangitis
5. Sepsis
6. Hyperamylasemia
7. Complications which are not currently known or observed may be present

CLINICAL BENEFIT

This device is used for the endoscopic removal of stones in the biliary system and foreign bodies, which can relieve biliary obstruction.

COMPATIBILITY

The device should be used with the endoscope with a working channel of at least 1.1mm

CAUTIONS

1. This device is intended for use by physicians or medical personnel who have received appropriate training, in accordance with the applicable laws. As a consumable of digestive endoscope, the product shall be used by professionals familiar with the operation technique of digestive endoscope.
2. This device is only allowed for clinical use in the endoscopy room.
3. Please read the instruction carefully before using it.

4. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

NOTES

A single copy of the instructions for use is provided in a case box (2 devices). If further copies are requested, please contact Micro-Tech, we will provide free of charge.

HOW SUPPLIED

The device is sterilized by EO and for single use only.

STORAGE

This device should be stored in a clean, well-ventilated, non-corrosive gas environment with a service life of 3 years.

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1. The device has been sterilized. Before use, carefully examine the instrument to make sure the product or aseptic package is not damaged during transportation. Do not use if there is any damage, please inform MT and send back the damaged products immediately.

WARNING: Please do not attempt to re-sterilize the device after use, which may result in infection control risk.

2. Confirm the validity of the product. Do not use if the product is beyond the expiry date.
3. Please refer to the applicable size of endoscope marked on the packing label to choose the appropriate product model.

Note: The device should be used with the endoscope with a working channel of at least 1.1mm.

4. Open the package and remove the device.
5. Operate the handle to check whether the basket can open and close smoothly. Do not use the basket if it cannot open and close smoothly.
6. Fully open the basket. Check the appearance of the device for breakage, excessive bending, fracture and distortion by visual inspection. Don't use if there is any damage.

WARNING: Please don't try to repair failed or damaged device.

WARNING: Please don't use device after disassembly.

【INSTRUCTIONS FOR USE】

Note: This device has no crushing stones function. In the event that the physician has decided that the calculus is big and it can't be extracted, please choose other device to crush the stones first.

1. After it is confirmed that the Extraction Basket is completely incorporated into the outer tube, insert the device into the channel of choledochoscope and slowly push it until the head end of device is exposed to the distal end of choledochoscope.
2. Advance the catheter beyond the stones or foreign bodies to be captured, slowly open the basket, grasp the stones or foreign bodies by pull the basket over them.

Note: Small movements and rotation of the device maybe necessary to grasp the targets.

3. Bring out the stones or foreign bodies from the bile duct by pull back the catheter.

Note: If the stones can't pass through the papilla, the endoscopic sphincterotomy is necessary.

When stone impaction, follow clinical and/or surgical standards of practice.

Note: Please don't retract retrieved stones and foreign bodies into the endoscope. This could damage the endoscope and/or device.

4. Pull the handle to close the basket.
5. Remove the device from the endoscope with the basket in the closed position.

PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

POST- PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Manufacturer		Apply to pancreaticobiliary system
	Date of manufacture		Use-by date
	Batch code		Catalogue number
	Not made with natural rubber latex.		Sterilized using ethylene oxide
	Do not re-sterilize		Do not use if package is damaged
	Keep away from sunlight		Keep dry
	Do not re-use		Consult instructions for use or consult electronic instructions for use
	Caution		Contents
	Open here		Sterile barrier system /sterile packaging
	Diameter		Compatible working channel
	Medical Device		Unique device identifier

【Packaging】 Packed in flexible peel pouch

【Production Date】 See label

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 3 Years

WARRANTY

Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. Any technical advice with respect to the product and guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp. GmbH (Europe)