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WARNINGS

1. **The product is intended for single use only!** DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed.
2. Do not use this instrument for any purpose other than its intended use.
3. The product is intended to be used in the adult and adolescent populations.

DEVICE NAME

Guidewire Locking Device

DEVICE DESCRIPTION

【CONTENTS】

Guidewire Locking Device is installed at the endoscope channel inlet. In the ERCP operation, the device and guide wire can be inserted into the endoscope through the sealing cap of the Guidewire Locking Device, up to two guide wires can pass the Guidewire Locking Device at the same time. When the device and guide wire enter the predetermined position, the guide wire can be locked on the Guidewire Locking Device, and then the exchange and operation of other working devices can be carried out.

【MODEL NUMBER】

The Guidewire Locking Device includes 2 specifications: MT-RGL-O and MT-RGL-P.

MT-RGL-O is compatible with Olympus/Fujinon endoscope; MT-RGL-P is compatible with Pentax endoscope.

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Table 1 Specifications of Guidewire Locking Device

No.	REF	Adaptive Endoscope
1	MT-RGL-O	Olympus/ Fujinon
2	MT-RGL-P	Pentax

【STRUCTURE】

The Guidewire Locking Device is composed with the guidewire locking assembly and the clamp seal assembly.


Figure 1 Product Drawing of Guidewire Locking Device

No.	Part name	No.	Part name
1	Left Lock Block	5	Seal Cap Top
2	Right Lock Block	6	Seal Cap Base
3	Antiskid Part	7	Spring
4	Shim	8	Slide Fixture Block

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by healthcare professional who received appropriate training in ERCP or with her/his supervision. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

INTENDED USE/INDICATIONS

The Guidewire Locking Device is an accessory to be used with endoscopic devices to lock the wire guide(s) in place during ERCP procedures.

CONTRAINDICATIONS

Including but not limited to: not suitable for ERCP procedures.

CAUTIONS

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
3. Please read the instructions for use entirely before use.
4. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

NOTES

1. The device is intended for single use only and single person only.
2. If the device cannot be used by sudden reason, please remove the device by right way.
3. A single copy of the instructions for use is provided in a case box (1 device). If further copies are requested, please contact Micro-Tech, we will provide free of charge.

HOW SUPPLIED

The package is supplied sterile by ethylene oxide gas. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

STORAGE

1. The product should be stored in a cool, dry, clean, well-ventilated, non-corrosive gas environment.
2. Do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

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1. Read the product label and select the appropriate device for the procedure and endoscope in use.
2. The contents are supplied STERILE. Inspect the pouch for damage, if damaged, do not use.
3. Open the package carefully after verifying the expiration date.
4. Carefully remove the device from its packaging.

Note: Do not use excessive force as this may damage the device and affect performance.

5. Visually inspect the device with particular attention to bends or breaks. If any component appears damaged, or not functioning correctly, do not use the device.

【INSTRUCTIONS FOR USE】

The Guidewire Locking Device should be fixed on the endoscope channel inlet prior to the ERCP procedure.

1. Fix the Guidewire Locking Device to the endoscope channel inlet. At the time of installation, press the slide fixture block of the Guidewire Locking Device, and then align the endoscope channel inlet, press the Guidewire Locking Device to the bottom, and then release the slide fixture block. Finally, make sure that the Guidewire Locking Device is securely fixed on the endoscope channel inlet (see Figure 2).


Figure 2


Figure 3

2. Insert the device through the sealing cap, align the device with the endoscope channel inlet and insert slowly in a straight manner. When the device and guide wire enter the predetermined position, the guide wire can be locked in the locking groove and the guide wire is locked (see Figure 3). If necessary, insert the end of the guide wire through the Antiskid Part to fix the guide wire.

3. Then, after the device is used, release the guide wire and exit the device.

PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Do not re-use		Do not resterilize
	Date of manufacture		Manufacturer
	Use-by date		Caution
	Catalogue number		Sterilized using ethylene oxide
	Keep away from sunlight		Authorized representative in the European Community/European Union
	Keep dry		Do not use if package is damaged and consult instructions for use
	Consult instructions for use or consult electronic instructions for use		Not made with natural rubber latex.
	Batch code		Sterile barrier system/ sterile packaging
	Medical device		Open here
	Importer		Unique device identifier

【Packaging】 Packed in flexible peel pouch

【Production Date】 See packaging

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Period of Validity】 3 Years

WARRANTY

1. **Limited Warranty to Buyer.** Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. In the event that Micro-Tech gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.
2. **Disclaimer and Release.** the warranties, obligations, and liabilities of micro-tech as set forth herein are exclusive. buyer hereby waives all other warranties, express or implied, arising by law or otherwise, with respect to the products and any other goods or services delivered by buyer, including, but not limited to: (1) all other express and implied warranties, including any implied warranty of merchantability or fitness for a particular purpose, and (2) any implied warranty arising from course of performance, course of dealing, or usage of trade.
3. **Implied Warranties.** The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.
4. **Limitation of Liability.** except to the extent prohibited by applicable law, micro-tech's liability under this warranty is limited to: (a) the replacement of the products or the re-supply of equivalent products; (b) the repair of the products or payment of the cost of repairing the products; or (c) payment of the cost of replacing the products or acquiring equivalent products. micro-tech shall have no obligation or liability, whether arising in contract (including warranty), tort (including active, passive, or imputed negligence, strict liability, or product liability) or otherwise, for any special, consequential, punitive, incidental, or indirect damages, or for loss of use, loss of revenue, loss of business, lost profit, or other financial loss arising out of or in connection with any product or other goods or services furnished by micro-tech, even if micro-tech was aware of the possibility of such damages or loss.

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Instruction page for electronic Instructions for Use (eIFU) website

To view / access IFUs, please go to <https://eifu.micro-tech.com.cn> by entering the product-specific IFU code RGL02.

The hardware and software requirements to view / access IFUs are as follows:

- Both the computer or mobile devices (mobile phones, pad, notebooks and computer) can be used, with current operating system (computer: Apple's macOS system, Microsoft's Windows system; Mobile: Apple iOS system, Google Android system, HarmonyOS system).
- The browser must be the latest generation of Google Chrome, Firefox, Safari or Microsoft Edge.

For a paper copy, provided in 7 days at no cost, please call + 49(0)211 73 27 626-0.

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Micro-Tech (Nanjing) Co., Ltd.
Micro-Tech Europe GmbH
Shanghai international Holding Corp. GmbH (Europe)