Intestinal Stent

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
Warnings

	The Intestinal Stent Delivery System is intended for single use only! DO NOT reuse, reprocess or re-sterilize the device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse or reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
	Use the stent system prior to the "Use By" date specified on the package.
	Pay attention to the instructions for use. Make sure to read the instructions before using.

Caution

Patients sensitive to Nickel Titanium (Nitinol) may suffer an allergic reaction to this implant. The stent should be used with caution and only after careful consideration in patients with significant preexisting pulmonary or cardiac disease. The device is intended for use by qualified endoscopists or radiology physicians who have received appropriate training. Radiographic equipment that provides high quality images is needed.

Informed consent should be obtained from all patients who undergo intestinal stent implant. The doctors must inform the patients of all the possible benefits and risks as well as the short term and long term complications related to the procedure. Because of the complexity of the diseases there may be other complications which are unpredictable or not listed that will lead to injury, illness or death of the patients.

Warranty

Micro-Tech Nanjing (MT) warrants that reasonable care and prudence has been exercised in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness.** Handling, storage, cleaning, and sterilization of this instrument as well as other factors related to the patient, diagnosis, treatment, surgical procedures, and other matters beyond MT's control directly affect the instrument and the results obtained from its use. MT's obligation under this warranty is limited to the repair or replacement of this instrument and MT shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MT neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **MT assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.**

Contents

Brief Introduction	6
I Device Name.....	6
II Description	6
III Indications for Use	7
IV Contraindications	7
V Potential Complications	7
VI Warnings	8
VII Precautions	8
Directions for Use	8
I Pre-procedure	8
Equipment prepare	8
Stent Preparation	8
Visual Inspection	8
Flush the Delivery System	9
Reconfirmation	9
II Procedure	9
Precaution	9
Direction for General procedure	10
Endoscopic Procedure	10
Interventional procedure.....	11
General procedure	12
III Post Procedure	13
Compatibility	13
How Supplied	14
Indications	14
Specifications	15
Contacts	15

Brief Introduction

I Device Name

Intestinal Stent

II Description

- The Intestinal Stent Delivery System comprises of two components: the implantable metallic stent and the delivery system (refer to Fig. 1a, 1b, 1c). The stent is made of Nitinol wire by weaving in a tubular mesh shape. This shape design can make the stent more flexible, compliant and self-expanding. The delivery system consists of three coaxial tubes. The outer tube serves to constrain the stent until being retracted during the stent deployment. Radiopaque marker bands situated on the inner and outer tubes aid in imaging during the deployment. The inner tube contains a central lumen that accommodates a 0.035 inch/0.89mm guide wire.

- There are two Intestinal Stent Delivery Systems here: Duodenal Stent Delivery System and Large Intestine Stent Delivery System, each of them will be elaborated respectively. And each system can be operated either OTW(Over the Wire) or TTS (Through the Scope). TTS Intestinal Stent Delivery System is compatible with minimum 4.2mm working channel of the endoscope.


Fig.1a OTW Intestinal Stent Delivery System


Fig.1b OTW Resealable Intestinal Stent Delivery System


Fig.1c TTS Intestinal Stent delivery system

Stent Characteristics

- The reason why the Nitinol is used as the material lies in its physical characteristics: excellent biocompatibility, prominent corrosion tolerance, shape memory effect and super elasticity. Initially, the stent is intenerated under the condition of 0~10 °C or in ice water where its shape can be changed, so the stent can be easily loaded to the delivery system. Under the conditions inside the human body, where the temperature is more than 33 °C, the stent will resume its original shape gradually after being deployed from the delivery system. The stent will engender a gentle radial force which acts on the inner wall of intestinal tract to expand the stricture gradually and rebuild the unobstructed lumen. Super elasticity under body temperature helps the stent accomodate the intestinal peristalsis. Because of the special designs of the stent, the patient will feel more comfortable by keeping the intestinal tract patent after the implant.

- Both ends of the stent are pliable and smooth enough without any sharp corners or burrs. This design can reduce risks of injury to the intestinal tract. (Refer to Fig. 2a)


Fig. 2a


Fig. 2b


Fig. 2c

- The stent can be covered by silicone on request. The silicone covering has characteristics of excellent biocompatibility and prominent corrosion tolerance. And it can inhibit the tumor ingrowth through the stent mesh as well as occlude the intestinal fistula. (Refer to Fig. 2b)

- Different kinds of stents are available for different situations. (Refer to Fig. 2c)

III Indications for Use

The Intestinal Stent is indicated for use in the treatment of intestinal strictures produced by malignant neoplasms, stricture in intestinal anastomotic stoma and the intestinal fistula occluding.

IV Contraindications

Contraindications Include, but not limited to:

- Patients with cardiac or pulmonary function failure
- Stricture caused by benign tumors
- Severe coagulopathy
- Severe gastrointestinal varices and bleeding caused by portal vein hypertension
- Severe intestinal adhesion and subsequent intestinal obstruction
- Bleeding caused by severe hemorrhoid and Perianal varices
- Bleeding caused by acute inflammation and cankerous colitis.

V Potential Complications

Complications related to the procedure include, but not limited to:

PROCEDURAL COMPLICATIONS

- Stent misplacement
- Intestinal perforation
- Infection

Bleeding

Pain

POST PROCEDURAL COMPLICATIONS

Stent occlusion due to tumor overgrowth at stent ends

Stent occlusion due to tumor ingrowth

Stent occlusion due to granulomatous tissue ingrowth

Restenosis due to granulomatous tissue formation at stent ends

Intestinal ulceration and/or perforation and/or hemorrhage

Stent Break

Stent migration

Recurrent of defecation difficulties related to stent occlusion or migration

Death

VI Warnings

Do not use if the pouch is opened or damaged before use.

Do not expose the delivery system to organic solvents (e.g. alcohol).

Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

VII Precautions

The delivery system is not designed for use with power injection systems.

Store in a cool, dry place.

Directions for Use

I Pre-procedure

Equipment prepare

Endoscope, flexible or rigid

A 0.035 inch (0.89mm) guide wire

Intestinal Stent Delivery System

Syringe for Irrigation

Dilation Balloon (as necessary)

Stent Preparation

- Select the right length of the stent:

Generally speaking, the stent should be 20~30mm longer than the stricture in length. The distal end of the stent should be 10mm-20mm below the lower margin of the stricture while the proximal end should be about 10 mm above the upper margin of the stricture.

- Select the right diameter of the stent:

Generally speaking, the diameter of the stent is around 18~30mm depending on the conditions of the strictures.

Notice: Do not release stent entirely, or it can't be reloaded. Please change another one in this situation, reloading without proper training could lead to operation failure and result in harm or danger to the patients.

Visual Inspection

- Open the outer package to inspect the pouch to make sure that it is free from the damage. Then carefully open the pouch and take the stent

delivery system from the tray. Make sure that the device is free from any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

Flush the Delivery System

- It is advisable to flush luer port or back handle (guide wire port) of the stent delivery system with a 10cc syringe of saline to expel air. (Refer to Fig. 3a, 3c) Continue to flush until the saline flows out of the distal catheter tip. (Refer to Fig. 3b, 3d).


Fig. 3a


Fig. 3b


Fig. 3c


Fig. 3d

Reconfirmation

- Observe the distal end of the outer tube to ensure that the stent is kept inside the outer tube. Do not use if the stent is partially deployed.

Notice:

- Before the procedure, carry radiography to define the location, diameter and length of the stricture so as to select a suitable stent.
- Patients should stop eating one day before the procedure. This device should be operated by qualified doctors in the procedure.
- Whether pre-dilation is necessary or not depends on the situation and physician's judgment. If the stent delivery system can pass the stricture with minimal difficulty, pre-dilation is not necessary. Perform pre-dilation only when the stricture is so tight enough that the stent delivery system can't traverse, because pre-dilation may increase the risks of perforation and migration.

II Procedure

Precaution

- There may be one gap between the proximal end of the stent and the delivery system because of transportation, which may result in difficulty in stent deployment. Some steps should be followed to eliminate the gap before being used in patients. First of all, insert one guide wire into the stent delivery system through the luer port and loose the safety lock. Then immobilize the front handle and push the back handle gently until the gap eliminates. Finally, lock the safety lock.

- Insert the delivery system slowly and carefully along the guide wire under some effectual monitoring.
- Patients treated by radiotherapy or chemotherapy may result in tumor shrinkage and subsequent stent migration. Radiotherapy or chemotherapy can be carried out 30 days later in order to greatly decrease the risk of stent migration.
- If the endoscope will be used to check the stricture, and the stricture is so serious that the endoscope cannot pass through, forcing the endoscope to pass without the radiographic guiding may cause perforation. It is better to use a stiff guide wire under X-ray, thus the possibility of perforation will be reduced greatly.
- Intraluminal circulating pressure and multiple corrosion by gastrointestinal secretion may lead to the metal fatigue and subsequent stent break.
- Patients should lie on their left side and the procedure should be monitored under endoscopic guidance with the aid of imaging equipments like X-RAY and CT.

Direction for General procedure

Endoscopic Procedure


Fig.4 Pass an endoscope into the duodenum.


Fig.5 Place a guide wire through the endoscope.


Fig.6 Perform pre-dilation if necessary.


Fig.7 Pass the endoscope through the stricture.


Fig.8 Hold the olive tip with one hand, pull the stylet out completely with the other hand.


Fig. 9 remove the endoscope and advance the intestinal stent delivery system over the guide wire into the stricture.


Fig.10a1 For Resheathable type: loose the safe lock first, then withdraw the front handle to visual marker to deploy part of the stent while immobilizing the back handle.

Fig.10a2 Withdraw the front handle to deploy the stent totally out of the delivery system.

Note: Before reaching the visual marker , the stent can be reloaded to the delivery system again. The process of pulling the stent back is that push forward the front handle.


Fig.10b For other types: Loose the safe lock first. Then withdraw the front handle by one hand to deploy the stent while immobilizing the back handle.


Fig.11 Remove the delivery system.


Fig.12 Confirm endoscopically that the stent has been completely deployed.

Interventional procedure


Fig.13 Place a floppy tip guide wire through the stricture.


Fig.14 Place an exchange catheter through the guide wire.


Fig.15 Exchange a floppy tip guide wire with a stiff guide wire (if necessary).


Fig.16 Perform pre-dilate if necessary.


Fig.17 Hold the olive tip with one hand, pull the stylus out completely with the other hand.


Fig.18 Advance the OTW intestinal Stent System over the guide wire towards the stricture.


Fig.19a1 For Resheathable type: loose the safe lock first, then withdraw the front handle to deploy part of the stent while immobilizing the back handle.

Fig.19a2 Withdraw the front handle to deploy the stent totally out of the delivery system.

Note: Before reaching the marker line , the stent can be reloaded to the delivery system again. The process of pulling the stent back is that push forward the front handle.


Fig.19b For other types: Loose the safe lock first. Then withdraw the front handle by one hand to deploy the stent while immobilizing the back handle.

Fig.20 Remove the delivery system.

Fig.21 Confirm radiographically that the stent has been completely deployed.

General procedure

1. Locate the stricture and perform pre-dilation if necessary.

Pass an endoscope through the mouth into the intestinal until the stricture. If necessary, perform pre-dilation by dilation balloon to a minimum diameter to accommodate the preloaded stent delivery catheter, or until the endoscope can pass. (Refer to Fig. 4 to 7).

2. Measure the stricture length and intestinal tract diameter.

If the endoscope will be used to measure the stricture, pass the endoscope through the stricture , the length of the stricture from the incisors and the stricture length are measured. And visually estimate the diameter of the esophageal tract around the stricture. More accurately, they can be measured by reviewing a recent CT scan around the stricture. Once the measurements have been confirmed, the size of the stent can be determined.

3. Insert the guide wire.

Place a 0.035 inch (0.89mm) guide wire through the endoscope and advance the guide wire along the intestinal tract into the distal end of the stricture in intestinal. Then remove the endoscope and make the appropriate stent ready.

4. Advance the OTW intestinal Stent System over the guide wire

pull the protective stylus out first. (Refer to Fig. 8)

Advance the OTW intestinal Stent System over the guide wire until the proximal end of the preloaded stent is above the upper margin of the stricture and the distal end of the stent is below the lower margin of the stricture. Confirm the position by the radiopaque markers on the stent with the aid of radiography.

5. Deploy the stent.

Loosen the safety lock first. Then retract the front handle by one hand to deploy the stent while immobilizing the back handle. (Refer to Fig.10)

6. Assess the deployed stent and remove the delivery system.

Confirm both endoscopically and radiographically that the stent has been completely deployed from the catheter and in the right position (Refer to Fig.12). Carefully remove the delivery system and the guide wire (Refer to Fig. 11).

III Post Procedure

After the procedure, carry out radiography to ascertain the position of the stent. Consult the doctors before taking any food. Follow-up examinations by radiography and endoscopy should be performed to check any signs of complications.

Compatibility

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that this stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely immediately after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial magnetic gradient of 720 Gauss/cm or less outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) ≤ 2.0 W/kg). ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).
- Non-clinical testing was conducted on the stent under the following conditions, and produced a maximum temperature rise of 3.3 °C
 - a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
 - a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.7 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR scanner.

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen of the Prostatic Stent or within approximately 5 mm of the position of the Prostatic Stent as found during non-clinical testing using the sequences: T1-weighted, spin echo pulse sequence and Gradient echo pulse sequence in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR system with body radiofrequency coil. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Additional Information: The safety of performing an MRI procedure in a patient with overlapping duodenal stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

How Supplied

The Intestinal Stent System is supplied sterile (by ethylene oxide gas) and is intended for SINGLE USE ONLY.

Indications

	Manufacturer		Authorized representative in the European community
	Do not reuse		Do not re-sterilize
	Keep away from sunlight		Keep dry
	Use by		Date of manufacture
	Catalogue number		Batch code
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Caution
	MR conditional		Working Length
	Diameter		Recommended Guide wire
	Open here		Sterile barrier system/sterile packaging
	Unique device identifier		

Specifications

【Name】 Intestinal Stent

【Components】 Titanium, Nickel and Silicone

【Structure】 Braided/Woven stent, covered with silicone covering

【Shape】 Cylinder with one or two ends are spherical or mushroom or cup or horn

【Scale】 Diameter: 18 to 30 mm; Length: 40 to 160mm

【Production Date】 See package

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Validity of the sterilization】 Two years

【Notes】 For single use only! If the pouch is damaged, the product can not be used

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MICRO-TECH (NANJING) CO., LTD.
Shanghai International Holding Corp. GmbH (Europa)
MICRO-TECH Europe GmbH