Biliary Stent

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
Warnings

	The Biliary Stent Delivery System is intended for single use only! DO NOT reuse, reprocess or re-sterilize the device. Reuse, reprocessing or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse or reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
	Use the stent system prior to the "Use By" date specified on the package.
	Pay attention to the instructions for use. Make sure to read the instructions before using.

Caution

Patients sensitive to Nickel Titanium (Nitinol) may suffer an allergic reaction to this implant. The stent should be used with caution and only after careful consideration in patients with significant preexisting pulmonary or cardiac disease. The device is intended for use by qualified endoscopists or radiology physicians who have received appropriate training. Radiographic equipment that provides high quality images is needed.

Informed consent should be obtained from all patients who undergo biliary stent implant. The doctors must inform the patients of all the possible benefits and risks as well as the short term and long term complications related to the procedure. Because of the complexity of the diseases there may be other complications which are unpredictable or not listed that will lead to injury, illness or death of the patients.

Warranty

Micro-Tech Nanjing (MT) warrants that reasonable care and prudence has been exercised in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness.** Handling, storage, cleaning, and sterilization of this instrument as well as other factors related to the patient, diagnosis, treatment, surgical procedures, and other matters beyond MT's control directly affect the instrument and the results obtained from its use. MT's obligation under this warranty is limited to the repair or replacement of this instrument and MT shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. MT neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **MT assumes no liability with respect to instruments reused, reprocessed or re-sterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for intended use, with respect to such instrument.**

Contents

Brief Introduction	6
I Device Name	6
II Description	6
III Indications for Use	7
IV Contraindications	7
V Potential Complications	7
VI Warnings	8
VII Precautions	8
Directions for Use	8
I Pre-procedure	8
Equipment prepare	8
Stent Preparation	8
Visual Inspection	8
Flush the Delivery System	9
Reconfirmation	9
II Procedure	9
Precaution	9
ERCP procedure	10
PTCD PROCEDURE	12
III Post Procedure	12
Compatibility	13
How Supplied	13
Indications	14
Specifications	14
Contacts	15

Brief Introduction

I Device Name

Biliary Stent

II Description

The Biliary Stent Delivery System comprises of two components: the implantable metallic stent and the delivery system (Refer to Fig. 1a,1b,1c and 1d). The stent is made of Nitinol wire by weaving in a tubular mesh shape. This shape design can make the stent more flexible, compliant and self-expanding. The delivery system consists of three coaxial tubes. The outer tube serves to constrain the stent until being retracted during the stent deployment. Radiopaque marker bands situated on the tubes and stents aid in imaging during the deployment. The interior tube contains a central lumen that accommodates a 0.035in./0.89mm guide wire.


Fig.1a TTS Biliary Stent Delivery System for ERCP


Fig.1b General Biliary Stent Delivery System for PTCD


Fig.1c TTS Resealable Biliary Stent Delivery System for ERCP


Fig.1d General Resheathable Biliary Stent Delivery System for PTCD

Stent Characteristics

The reason why the Nitinol is used as the material lies in its physical characteristics: excellent biocompatibility, prominent corrosion tolerance, shape memory effect and super elasticity. Initially, the stent is intenerated under the condition of 0~10 °C or in ice water where its shape can be changed, so the stent can be easily loaded to the delivery system. Under the conditions inside the human body, where the temperature is more than 33 °C, the stent will resume its original shape gradually after being deployed from the delivery system. The stent will engender a gentle radial force which acts on the inner wall of biliary duct to expand the stricture gradually and rebuild the unobstructed lumen. Because of the special designs of the stent, the patient will feel more comfortable by keeping the biliary duct patent after the implant.

Both ends of the stent are pliable and smooth without any sharp corners or burrs. This design can reduce risks of injury to the biliary duct.(Refer to Fig. 2a)


Fig. 2a


Fig. 2b

The stent can be covered by silicone on request. The silicone covering has characteristics of excellent biocompatibility and prominent corrosion tolerance to gastric acid. And it can inhibit the tumor ingrowth through the stent mesh as well as occlude the biliary duct fistula. (Refer to Fig. 2b)

III Indications for Use

The Biliary Stent is indicated for use in the treatment of biliary strictures caused by malignant neoplasms.

IV Contraindications

Contraindications include, but not limited to:

Patients with obvious tendency of haemorrhage

Patients allergic to iodine and the anesthetic drugs;

Patients with hydroperitoneum and liver separation from the abdominal wall;

Patients with space occupying lesion in puncturing path;

V Potential Complications

Complications related to the procedure include, but not limited to:

PROCEDURAL COMPLICATIONS

Stent misplacement
Perforation
Infection
Bleeding
Pain

POST PROCEDURAL COMPLICATIONS

Stent occlusion due to sludge
Ulceration
Perforation
Hemorrhage
Stent Break
Stent migration
Stent occlusion due to tumor ingrowth
Stent occlusion due to granulosomatous tissue ingrowth
Recurrent obstructive dysphagia due to stent occlusion or migration
Death

VI Warnings

Do not use if the pouch is opened or damaged before use.
Do not expose the delivery system to organic solvents (e.g. alcohol).
Upon completion of procedure, dispose of device per institutional guidelines for biohazardous medical waste.

VII Precautions

The delivery system is not designed for use with power injection systems.
Store in a cool, dry place.

Directions for Use

I Pre-procedure

Equipment Required

Endoscope, flexible
A 0.035" (0.89mm) Guide Wire
Biliary Stent Delivery System
Syringe for irrigation
Dilatation Balloon (as necessary)

Stent Preparation

- Select the length of the stent:

Generally speaking, the stent should be 20mm longer than the stricture in length. The distal end of the stent should be 10 mm below the lower margin of the stricture while the proximal end should be about 10 mm above the upper margin of the stricture.

- Select the diameter of the stent:

Generally speaking, the diameter of the stent is around 8~10mm depending on the conditions of the strictures. 6-8mm is advisable if two biliary stents need to be implanted.

Visual Inspection

- Open the outer package to inspect the pouch to make sure that it is free

from the damage. Then carefully open the pouch and take the stent delivery system from the tray. Make sure that the device is free from any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

Flush the Delivery System

- It is advisable to flush back handle (guide wire port) of the stent delivery system with a 10cc syringe of saline to expel air. (Refer to Fig. 3a) Continue to flush until the saline flows out of the distal catheter tip. (Refer to Fig. 3b)


Fig. 3a


Fig. 3b

Reconfirmation

- Observe the distal end of the outer tube to ensure that the stent is totally placed inside the outer tube. Do not use if the stent is partially deployed.

II Procedure

Precaution

- Insert the delivery system slowly and carefully along the guide wire under some effectual monitoring.
- Patients treated by radiotherapy or chemotherapy may result in tumor shrinkage and subsequent stent migration. Radiotherapy or chemotherapy can be carried out 30 days later in order to greatly decrease the risk of stent migration.
- Nosogenetic complexities and intraluminal circulating pressure may lead to the metal fatigue and subsequent stent break.
- Sphincterotomy is not always essential for stent delivery, but may be performed based on the physicians' judgments.
- Biliary stent can be performed by ERCP or PTCD. Both of them will be elaborated respectively. ERCP Biliary Stent Delivery System is compatible with minimum 3.2mm working channel for a 7.5Fr delivery system or minimum 3.7 mm working channel for a 8.5Fr delivery system.

ERCP procedure

1. Define the stricture and biliary duct diameter.

Pass an endoscope until the duodenal papilla. Then advance the cannulating catheter through the working channel of the endoscope. Perform the catheter to pass through the papilla and inject contrast agent. Observe the contrast agent by radiograph to define the stricture and biliary duct diameter.

2. Insert the guide wire and perform pre-dilation if necessary.

Place a 0.035 inch (0.89mm) guide wire through the catheter until the biliary duct, then remove the catheter. Perform sphincterotomy if necessary. If the stricture is so tight, perform pre-dilation by dilation balloon to accommodate the stent delivery system.

3. Advance the Biliary Stent Delivery System over the guide wire

Pull protective stylus out first.

Advance The Biliary Stent Delivery System over the guide wire into the biliary

duct until the proximal end of the preloaded stent is above the upper margin of the stricture and the distal end of the stent is below the lower margin of the stricture. Confirm the position by the radiopaque markers on the stent with the aid of radiography.

4. Deploy the stent.

General procedure: Loose the safe lock first. Then retract the front handle by one hand to deploy the stent while immobilizing the back handle by the other hand.

5. Assess the deployed stent and remove the delivery system.

Confirm radiographically that the stent has been completely deployed from the system and in the right position. Carefully remove the delivery system and the guide wire while watching not to dislodge the stent with the olive tip.

ERCP Procedure


Fig.4 Pass an endoscope until the duodenal papilla.

Fig.5 Calculate through the papilla.


Fig.6 Inject contrast agent to define the stricture and biliary duct diameter.


Fig.7 place a 0.035 inch (0.89mm) guide wire through Endoscope with the visualization of radiography.

Fig.8 Perform sphincterotomy if necessary.


Fig.9 Perform pre-dilation if necessary.


Fig.10 Hold the olive tip with one hand, pull the stylus out completely with the other hand.


Fig.11 Estimate the stricture length by measuring length of the guide wire inserted when withdraw it.


Fig.12 Advance the Biliary Stent Delivery System over the guide wire


Fig.13 Loose the safe lock first. Then retract the front handle by one hand to deploy the stent while immobilizing the back handle by the other hand.


Fig.14a For Resheathable type: loose the safe lock first, then withdraw the front handle to visual marker to deploy part of the stent while immobilizing the back handle.

Fig.14b Withdraw the front handle to deploy the stent totally out of the delivery system.

Note: Before reaching the visual marker , the stent can be reloaded to the delivery system again. The process of pulling the stent back is that push forward the front handle.


Fig.15 Remove the delivery system.


Fig.16 Confirm radiographically that the stent has been completely deployed.

PTCD PROCEDURE

1. Take the operation of PTCD by puncturing through the left or right biliary duct. After the guide wire and catheter passed the stricture and settled in the duodenum, insert the sheath of the catheter. Take the radiography examination to confirm the length of the stenosis and select a suitable stent. The selected stent should be 20mm longer than the stenosis.
2. Perform pre-dilate if necessary.
3. Under the X-ray, insert the delivery system into the biliary tract along the guide wire till the stricture. When the proximal end of the stent is about 5~10mm higher than the proximal margin of the stent, deploy the stent by retracting the front handle by one hand while immobilizing the back handle by the other hand.
4. Confirm radiographically that the stent has been completely deployed from the system and in the right position.
Use a balloon to help the stent expand if necessary.
5. Carefully remove the delivery system and the guide wire while watching not to dislodge the stent with the olive tip.

PTCD Procedure


Fig.17 Perform PTCD by puncturing into the left or right biliary branch.


Fig.18 Confirm the puncture by using a syringe to collect a small amount of bile.


Fig.19 Radiographically define the stricture by injecting contrast media through the catheter


Fig.20 Insert the guide wire to the distal end of the stricture.


Fig.21 Perform pre-dilation if necessary.


Fig.22 Measure the stricture length while pulling back the guide wire or inspecting the scale on the sheath.


Fig.23 Advance the delivery system over the guide wire


Fig.24 Deploy the stent


Fig.25 Remove the delivery system.


Fig.26 Confirm radiographically that the stent is completely deployed.

III Post Procedure

If the stricture reoccurs after stent implantation, it is unfavorable to dilate the stricture by balloon.

Compatibility

This symbol means the device is MR Conditional.

Non-clinical testing has demonstrated that this stent is MR Conditional according to ASTM F2503. A patient with this stent can be scanned safely immediately after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3 Tesla or less
- Maximum spatial magnetic gradient of 720 Gauss/cm or less outside of scanner covering, accessible to a patient or individual.

MRI-Related Heating

- 1.5 and 3.0 Tesla Systems: It is recommended to scan in normal operation mode (whole body averaged specific absorption rate (SAR) $\leq 2.0\text{W/kg}$). ("Normal Operating Mode" is defined as the mode of operation of the MR system in which none of the outputs have a value that cause physiological stress to the patient) for 15 minutes of scanning (i.e., per scanning sequence).
- Non-clinical testing was conducted on the stent under the following conditions, and produced a maximum temperature rise of 2.8 °C
 - a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.1 W/kg) for 15 minutes of MR scanning in a 1.5 Tesla Magnetom (Siemens Medical Solutions, Malvern, PA, Software Numaris/4) MR scanner.
 - a maximum whole body average specific absorption rate (SAR) of 2.9 W/kg (corresponding to a calorimetry measured value of 2.7 W/kg) for 15 minutes of MR scanning in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR scanner.

Image Artifacts

MR image quality may be compromised if the area of interest is within the lumen of the Biliary Stent or within approximately 5 mm of the position of the Biliary Stent as found during non-clinical testing using the sequences: T1-weighted, spin echo pulse sequence and Gradient echo pulse sequence in a 3.0 Tesla Excite (GE Electric Healthcare, Milwaukee, WI, Software 14X.M5) MR system with body radiofrequency coil. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

Additional Information: The safety of performing an MRI procedure in a patient with overlapping duodenal stents or other MRI-conditional device(s) in direct contact with this device has not been determined. Performing MRI in such situations is not recommended.

How Supplied

The Biliary Stent System is supplied sterilized (by ethylene oxide) and is intended for SINGLE USE ONLY.

Indications

	Manufacturer		Authorized representative in the European community
	Do not reuse		Do not re-sterilize
	Keep away from sunlight		Keep dry
	Use by		Date of manufacture
	Catalogue number		Batch code
	Consult instructions for use or consult electronic instructions for use		Do not use if package is damaged
	Sterilized using ethylene oxide		Caution
	MR conditional		Working Length
	Diameter		Recommended Guide wire
	Open here		Sterile barrier system/sterile packaging
	Unique device identifier		

Specifications

【Name】 Biliary Stent

【Components】 Titanium, Nickel and Silicone

【Structure】 Braided/Woven stent, covered with silicone covering

【Shape】 Cylinder with one or two horn ends

【Scale】 Diameter: 6 to 10 mm; Length: 40 to 120mm

【Production Date】 See package

【Sterilization】 Sterilized by EO (ethylene oxide) gas

【Validity of the sterilization】 Two years

【Notes】 For single use only! If the pouch is damaged, the product can not be used

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MICRO-TECH (NANJING) CO., LTD.
Shanghai International Holding Corp. GmbH (Europa)
MICRO-TECH Europe GmbH