Multiple Band Ligator Set Instructions for Use
Mehrband-Ligator-Set Gebrauchsanweisung für
Ensemble de ligateur multibandes Mode d'emploi
Conjunto de Ligadores de Multibanda Instrucciones de uso
Set di Legatura Multibanda Istruzioni per l'uso
Σετ Απολίνωσης Πολλαπλών Δακτυλίων Οδηγίες Χρήσης
Flerbånds ligatorsæt Brugsanvisning
Conjunto de Ligador Multibanda Instruções de utilização
Çoklu Bantlı Ligator Seti Kullanım Talimatları
Комплект с многолентови лигатори Инструкции за употреба
Set vícekroužkových ligátorů Návod k použití
Moninauhainen ligaattorisarja käyttöohjeet
Többcsatornás ligátorszett Használati útmutató
Kelių juostų ligatorių rinkinys naudojimo instrukcija
Vairāku joslu ligatoru komplekts Lietošanas instrukcija
Flerbånds ligator sett Bruksanvisning
Set ligator cu bandă multiplă Instrucțiuni de utilizare
Viacpásmová ligátorová súprava Návod na použitie
Komplet ligatorja z več trakovi Navodila za uporabo
Uppsättning av bandligering Bruksanvisning
Zestaw ligatorów wielopasmowych Instrukcje obsługi
Ligator Set met meerdere banden Gebruiksaanwijzing van
Multiple Band Ligator komplekt Kasutusjuhend
Komplet višepojasnih ligatora Upute za uporabu

Contents

Multiple Band Ligator Set Instructions for Use	1
Mehrband-Ligator-Set Gebrauchsanweisung für	8
Ensemble de ligateur multibandes Mode d'emploi	16
Conjunto de ligadores de multibanda Instrucciones de uso	24
Set di Legatura multibanda Istruzioni per l'uso	32
ΣΕΤ Απολίνωσης Πολλαπλών Δακτυλίων Οδηγίες Χρήσης	40
Flerbånds ligatorsæt Brugsanvisning	48
Conjunto de Ligador Multibanda Instruções de utilização	56
Çoklu Bantlı Ligator Seti Kullanım Talimatları	64
КОМПЛЕКТ С МНОГОЛЕНТОВИ ЛИГАТОРИ Инструкции за употреба	72
Set vícekroužkových ligátorů Návod k použití	80
Moninauhainen ligaattorisarja käyttöohjeet	88
Többszatosornás ligátorszett Használati útmutató	96
Kelių juostų ligatorių rinkinys naudojimo instrukcija	104
vairāku joslu ligatoru komplekts Lietošanas instrukcija	112
Flerbånds ligator sett Bruksanvisning	120
Set ligator cu bandă multiplă Instrucțiuni de utilizare	127
Viacpásmová ligátorová súprava Návod na použitie	135
Komplet ligatorja z več trakovi Navodila za uporabo	143
Uppsättning av bandligering Bruksanvisning	150
Zestaw ligatorów wielopasmowych Instrukcje obsługi	157
Ligator Set met meerdere banden Gebruiksaanwijzing van	165
Multiple Band Ligator komplet Kasutusjuhend	173
Komplet višepojasnih ligatora Upute za uporabu	181
Instruction page for Electronic instructions for use website	189

WARNING

1. **The product is intended for single use only!** DO NOT re-use, re-sterilize, and/or reprocess. Re-use, re-sterilization or reprocessing may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Re-use, re-sterilization or reprocessing may also create a risk of contamination of the device and/or cause patient infectious disease(s). Contamination of the device may lead to injury, illness or death of the patient. Micro-Tech assumes no liability with respect to instruments reused, re-sterilized or reprocessed.
2. Do not use this instrument for any purpose other than its intended use.
3. Do not insert the instrument into the endoscope unless there is a clear endoscopic field of view. If the distal end of the insertion portion in the endoscopic field of view cannot be seen, do not use the device. Insertion without clear endoscopic field of view could cause patient injury, such as perforation, hemorrhage or mucous membrane damage. It may also damage the endoscope and/or instrument.
4. The product is intended to be used in the adult and adolescent populations.
5. This device is available in two types: one is not made with natural rubber latex, the other is made with natural rubber latex. Please read the labeling carefully before use to determine the appropriate device for use.
6. Use of latex ligation bands may be contraindicated in patients with known hypersensitivity to latex.

DEVICE NAME

【COMMON NAME】

Multiple Band Ligator Set

DEVICE DESCRIPTION

【CATALOGUE NUMBER】

MBLS-X-Y Z-NL-T

	Set type: blank: Set, T: Non-Set
	Band type: blank: latex band, NL: Non-latex band
	Market Code: blank: China, F: Countries except China
	Bands Quantity
	OD of Endoscope: Blank (9.4-13.0)mm
	XL-Extra Large (11.0-14.0)mm
	Product Code: MBLS-Multiple Band Ligator Set

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The specifications of the product are represented by Endoscope OD, Band Quantity, Band Type and Set type, as indicated in Table 1.

Table 1 Specification and parameter


Endoscope OD ϕ (mm)	Band Quantity (pcs)	Band Type	Set type
9.4~14	3, 4, 6, 7	Non-Latex, Latex	Set, Non-Set

【STRUCTURE】

The Multiple band Ligator Set is mainly composed of set and non-set. The set is composed of Ligator Cartridge with preloaded bands, Handle Unit, Pull Catheter Set and Irrigation Needle, while the non-set is composed of Ligator Cartridge with preloaded bands and Pull Catheter Set.

1. The Ligator Cartridge with preloaded bands ---- This is a transparent cylinder with an elastic Soft Cap that allows it to be fixed on the distal end of the endoscope. The bands are mounted on Barrel in a stretched condition, each band is blue (yellow) in color with the second from last band being neutral in color. A Trigger Thread with Trigger Beads is placed under the bands. This Trigger Thread also has a knot at the proximal end to attach to the Pull Thread. The Pull Thread is passed through the working channel of the endoscope and attaches


to the Handle Unit. The bands are fashioned to be released by pulling the Pull Thread using the handle.

1. Trigger Thread; 2. Soft Cap; 3. Band; 4. Barrel; 5. Trigger Bead; 6. Pull Thread

Pic 1. The Ligator Cartridge with preloaded bands


2. The Handle Unit ---- When assembled, rotating the Knob of the Handle Unit, in a clockwise direction, pulls the Pull Thread and wraps it around the Winding Drum. Pulling the Pull Thread results in band deployment. To assist in set up, pulling the Knob away from the Handle Unit, permits free clockwise or anti-clockwise rotation. Pushing the Knob towards the Handle Unit, locks the rotation for use.

7. Knob; 8. Handle Base; 9. Winding Drum; 10. Mounting Base;

Pic 2. The Handle Unit

3. Pull Catheter Set ---- This Pull Catheter has a Hook on either end. It is used to pull the Pull Thread of the Ligator Cartridge through the working channel of the endoscope for attachment to the Handle Unit.

11. Hook; 12. Pull Catheter

Pic 3. Pull Catheter Set

4. Irrigation Needle ---- This may be attached to an irrigation syringe. The needle can be inserted through the Valve Opening in the handle to provide irrigation as required.

Pic 4. Irrigation Needle

【USER INFORMATION/TRAINING/QUALIFICATIONS】

This device is intended for use by healthcare professional who received appropriate training in digestive endoscopy or with her/his supervision. As an accessory of digestive endoscopy, the product shall be used by professionals familiar with the operation technique of digestive endoscopy.

INTENDED USE/ INDICATIONS

The Multiple Band Ligator Set is used to endoscopically ligate esophageal varices at or above the gastroesophageal junction or to ligate internal hemorrhoids.

CLINICAL BENEFIT STATEMENT

The device can achieve prevention of bleeding and effective hemostasis, while effectively eliminating esophageal varices and the symptoms of internal hemorrhoids.

PERFORMANCE CHARACTERISTICS

Multiple Band Ligator Set is suitable for endoscopes of 9.4 to 14.0 mm in diameter.

Multiple Band Ligator Set has a ligation-design with or without latex.

CONTRAINDICATIONS

For esophageal varices:

1.Contraindications to upper gastrointestinal endoscopy, 2.Severe cardiopulmonary disease, 3.Severe coagulopathy, 4.Severe liver and kidney disease, massive ascites, 5.Severe hepatic encephalopathy, 6.Uncorrected hemorrhagic shock, 7.Use of latex ligation bands may be contraindicated in patients with known hypersensitivity to latex, 8. Any other conditions the physician judges unsuitable for use.

For internal hemorrhoids:

1.Severe cardiopulmonary disease, 2.Severe liver and kidney disease, 3.Severe coagulopathy, 4.Severe inflammatory bowel disease, 5.Rectal malignant tumor, 6.Anal fissure, perianal abscess, 7.Pregnancy or menstrual period, 8.Use of latex ligation bands may be contraindicated in patients with known hypersensitivity to latex, 9.Any other conditions the physician judges unsuitable for use.

POTENTIAL COMPLICATIONS

For esophageal varices:

1.Retrosternal pain, 2.Stricture, 3.Ulceration, 4.Dysphagia, 5.Hemorrhage, 6.Infection, 7.Perforation, 8.Bands slipping off. Complications which are not currently known or observed may be present.

For internal hemorrhoids:

1.Pain, 2.Hemorrhage, 3.Urinary symptoms, 4.Edema, 5.Bulge, 6.Ulceration, 7.Infection, 8.Bands slipping off, 9.Anal incontinence. Complications which are not currently known or observed may be present.

CAUTIONS

1. A thorough understanding of the technical principles, clinical applications, and associated risks is expected before usage.
2. Inform the patient of the details of the operation and all of the potential risks and complications, which may lead to injury, illness, or death, and ensure that the patient has expressed his/her acceptance.
3. Please read the instructions for use in their entirety before use.

STORAGE

1. Storage conditions: in a room with a temperature range of 10-25℃ and a humidity range of $\leq 80\%$.
2. Do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation.

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
1. Read the product label and select the appropriate device for the procedure and endoscope in use.
2. The contents are supplied STERILE. Inspect the pouch for damage, do not use if the package is damaged.
3. Open the package carefully after verifying the expiration date.
4. Carefully remove the device from its packaging and uncoil it.

Note: Do not use excessive force as this may damage the device and affect performance.

5. Visually inspect the device with particular attention to bends or breaks. If it appears damaged, do not use

the device.

6. Feed one of the hooked ends of the Pull Catheter through the Mounting Base of the Handle Unit. Push the catheter until it emerges through the Valve Opening.

7. Assemble the endoscope and open the endoscopes biopsy valve cover. Insert the other end of the Pull Catheter into the biopsy channel of the endoscope until it emerges out of the distal end.
 8. Insert the Mounting Base of the Handle Unit directly through the biopsy valve into the instrument channel. Ensure that the handle is firmly seated and pull the Knob away from the Handle Unit into the free rotation position.
 9. Attach the free end of the Pull Thread from the Ligator Cartridge to the hook at distal end of the Pull Catheter. Ensure that the knot at the end of the Pull Thread is securely ensnared in the hook.
 10. Gently withdraw the Pull Catheter from the handle end and draw the Trigger Thread through the Valve Opening.
 11. Slip the proximal end of the Pull Thread into the groove on the Winding Drum of the Handle. Pull gently on the Pull Thread so that the knot disappears into the Tie Hole and the Pull Thread is firmly anchored to the Winding Drum.
 12. Push the knob in towards the Handle Unit to engage the rotation lock. Start rotating the Knob in a clockwise direction to take up the slack in the thread while placing the Ligator Cartridge close to the distal end of the scope. Stop rotating when the cartridge touches the scope.
 13. Secure the cartridge onto the distal end of the endoscope by pushing the Soft Cap over the tip of the endoscope.
- Note:** Ensure the cartridge is securely attached to the endoscope prior to proceeding.
14. Rotate the Knob to take up the rest of the slack, until the distal tip of the endoscope begins to move slightly. This indicates that all the slack has been taken up and the device is now ready for use.
 15. Attach a water-filled syringe to the Irrigation Needle, to be used as required.

【INSTRUCTIONS FOR USE】

For esophageal varices:

1. Introduce the endoscope, with Ligator Cartridge attached, into the alimentary canal.
2. Visualize the target esophageal varices.
3. Press the distal tip of the Ligator Cartridge against the esophageal varices.
4. Once complete contact is made around the esophageal varices, activate the suction of the endoscope and draw the tissue into the barrel of the Ligator Cartridge. When the tissue is completely sucked into the barrel, a “red out” will be observed.
5. Slowly rotate the Knob on the Handle Unit clockwise until a band is released. The tension on the Pull Thread / Knob of the Handle Unit will reduce immediately after the band has been released and deployed.

Note: Usually a band will deploy by turning the Knob a little more than half a turn. If a band does not deploy even after a full turn, pull the Knob away from the Handle Unit to release the rotation lock. Ease the tension on the Pull Thread and then push the Knob towards the Handle Unit to the locked position. Repeat step 5

6. Release suction and confirm band placement in the endoscopic view.
7. Rotate the Knob of the Handle Unit clockwise to take up any residual slack in the Pull Thread. To ligate

another esophageal varices, visualize a new target and repeat steps 3, 4, 5, 6 and 7.

Note: A neutral color band indicates that one more ligation band is available for use.


For internal hemorrhoids:

1. Introduce the endoscope, with Ligator Cartridge attached, into the rectum.
2. Visualize the target internal hemorrhoid.
3. Press the distal tip of the Ligator Cartridge against the internal hemorrhoid.
4. Once complete contact is made around the internal hemorrhoid, activate the suction of the endoscope and draw the tissue into the barrel of the Ligator Cartridge. When the tissue is completely sucked into the barrel, a “red out” will be observed.
5. Slowly rotate the Knob on the Handle Unit clockwise until a band is released. The tension on the Pull Thread / Knob of the Handle Unit will reduce immediately after the band has been released and deployed.

Note: Usually a band will deploy by turning the Knob a little more than half a turn. If a band does not deploy even after a full turn, pull the Knob away from the Handle Unit to release the rotation lock. Ease the tension on the Pull Thread and then push the Knob towards the Handle Unit to the locked position. Repeat step 5

6. Release suction and confirm band placement in the endoscopic view.
7. Rotate the Knob of the Handle Unit clockwise to take up any residual slack in the Pull Thread. To ligate another hemorrhoid, visualize a new target and repeat steps 3, 4, 5, 6 and 7.

Note: A neutral color band indicates that one more ligation band is available for use.

PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative, and/or local government policy.

POST-PROCEDURE

Any serious incident that occurs in relation to this device should be reported to the manufacturer and relevant local regulatory authority.

SYMBOL INDICATIONS

	Do not re-use		Do not re-sterilize
	Date of Manufacture		Manufacturer
	Use-by date		Caution
	Contains or presence of natural rubber latex		Not made with natural rubber latex.
	Catalogue Number		Sterilized using ethylene oxide
	Authorized representative in the European Community/ European Union		Keep away from sunlight
	Keep dry		Do not use if package is damaged and consult instructions for use
	Consult instructions for use or consult electronic instructions for use		Compatible working channel
	Diameter		Working length
	Batch code		Contents
	Sterile barrier system/sterile packaging		Medical Device
	Open here		Open here
	Importer		Unique device identifier
	Ambient Temperature		Relative humidity

【 **Packaging** 】 Preformed rigid tray with a die-cut lid or flexible peel pouch.

【 **Production Date** 】 See packaging

【 **Sterilization** 】 Sterilized by EO (ethylene oxide) gas

【 **Period of Validity** 】 1 year

【 **Compatible Working Channel** 】 $\geq \varnothing 2.8$ mm

【 **Endoscope Diameter** 】 9.4 ~ 14.0 mm

【 **Applicable Endoscopes** 】 The device is to be used with CE Marked endoscopes.

WARRANTY

1. Limited Warranty to Buyer. Micro-Tech warrants to Buyer that, for the earlier of one (1) year from the date of purchase, or until the product is used by Buyer, the products will be free from defects in materials and workmanship when stored and used in accordance with the instructions for storage and use provided by Micro-Tech and in accordance with applicable regulatory requirements. Descriptions or specifications appearing in Micro-Tech's literature are meant to generally describe the products and do not constitute any express warranties. In the event that Micro-Tech gives technical advice with respect to the product, it is agreed that such advice is given without any liability on Micro-Tech's part. Any guarantee of specific properties of or in the products shall only be effective if and to the extent specifically confirmed by Micro-Tech in writing. These warranties shall not apply for product failure or deficiency due to improper storage, alteration, or the consequences of uses for which the products were not designed or that adversely affect the products' integrity, reliability, or performance.
2. Disclaimer and Release. **THE WARRANTIES, OBLIGATIONS, AND LIABILITIES OF Micro-Tech AS SET FORTH HEREIN ARE EXCLUSIVE. BUYER HEREBY WAIVES ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, ARISING BY LAW OR OTHERWISE, WITH RESPECT TO THE PRODUCTS AND ANY OTHER GOODS OR SERVICES DELIVERED BY BUYER, INCLUDING, BUT NOT LIMITED TO: (1) ALL OTHER EXPRESS AND IMPLIED WARRANTIES, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, AND (2) ANY IMPLIED WARRANTY ARISING FROM COURSE OF PERFORMANCE, COURSE OF DEALING, OR USAGE OF TRADE.**
3. Implied Warranties. The purchase of products may be subject to laws in the territories applicable to the sale of the products by Micro-Tech to Buyer that impose implied warranties, conditions, or obligations that cannot be excluded, restricted, or modified, or can only be excluded, restricted, or modified to a limited extent. The provisions of Paragraphs 2 and 4 shall apply to the greatest extent allowed by such laws.
4. Limitation of Liability. **EXCEPT TO THE EXTENT PROHIBITED BY APPLICABLE LAW, Micro-Tech'S LIABILITY UNDER THIS WARRANTY IS LIMITED TO: (a) THE REPLACEMENT OF THE PRODUCTS OR THE RE-SUPPLY OF EQUIVALENT PRODUCTS; (b) THE REPAIR OF THE PRODUCTS OR PAYMENT OF THE COST OF REPAIRING THE PRODUCTS; or (c) PAYMENT OF THE COST OF REPLACING THE PRODUCTS OR ACQUIRING EQUIVALENT PRODUCTS. Micro-Tech SHALL HAVE NO OBLIGATION OR LIABILITY, WHETHER ARISING IN CONTRACT (INCLUDING WARRANTY), TORT (INCLUDING ACTIVE, PASSIVE, OR IMPUTED NEGLIGENCE, STRICT LIABILITY, OR PRODUCT LIABILITY) OR OTHERWISE, FOR ANY SPECIAL, CONSEQUENTIAL, PUNITIVE, INCIDENTAL, OR INDIRECT DAMAGES, OR FOR LOSS OF USE, LOSS OF REVENUE, LOSS OF BUSINESS, LOST PROFIT, OR OTHER FINANCIAL LOSS ARISING OUT OF OR IN CONNECTION WITH ANY PRODUCT OR OTHER GOODS OR SERVICES FURNISHED BY Micro-Tech, EVEN IF Micro-Tech WAS AWARE OF THE POSSIBILITY OF SUCH DAMAGES OR LOSS.**

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Micro-Tech (Nanjing) Co., Ltd.
MICRO-TECH Europe GmbH
Shanghai International Holding Corp. GmbH(Europe)